
ETHICON Proximate™

en Reloadable Vascular Linear Stapler and Reload

- cs Opakovaně naplnitelný vaskulární lineární stapler a náplně
da Vaskulær lineær stapler med genpåfyldningsfunktion og magasin
de Nachladbarer vaskulärer Linear Stapler und Magazin
el Επαναφορτώσιμο ευθύγραμμο συρραπτικό για αγγεία και ανταλλακτική κασέτα
es Grapadora lineal vascular recargable y recarga
fi Uudelleenladattava vaskulaarinen suorasulkuinstrumentti ja vaihtokasetti
fr Agrafeuse linéaire vasculaire rechargeable et chargeur
hu Újratölthető érendszeri lineáris kapcsozóműszer és utántöltő
it Suturatrice lineare vascolare ricaricabile e ricarica
nl Herlaadbare vasculaire lineaire stapler en vulling
pl Liniowy stapler naczyniowy do ponownego ładowania i magazynek
pt Agrafador linear vascular recarregável e recargas
sv Omladdningsbar linjär vaskulär stapler och magasin
tr Tekrar Yüklenebilir Vasküler Lineer Stapler ve Yedek Ünite
zh 可重复装钉直线型血管吻合器和钉仓
bg, et, hr, id, lt, lv, mk, no, ro, ru, sk, sl, sr, vi: www.e-ifu.com










Please read all information carefully.

Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.

Important: This package insert is designed to provide instructions for use of the PROXIMATE™ Reloadable Vascular Linear Stapler and Reload. It is not a reference to surgical stapling techniques.

ETHICON™

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Indications

The PROXIMATE™ Reloadable Vascular Linear Stapler has application in internal tissue which can easily be compressed to 1 mm in thickness and may also be used to ligate pulmonary vessels.

Contraindications

- Do not use the instrument with white reload on any tissue that requires excessive force to compress to 1.0 mm or on any tissue that compresses easily to below 1.0 mm.
- Do not use the instrument on ischemic or necrotic tissue.
- Do not use the instrument on the aorta.
- Do not use the instrument on solid organs, such as liver or spleen, where attempted compression would be destructive.
- This instrument is not intended for use when surgical stapling is contraindicated.

Undesirable Side Effects / Residual Risks

Undesirable side effects and risks associated with surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures or damaged devices.

Device Description / Performance Characteristics

The PROXIMATE™ Reloadable Vascular Linear Stapler and Reload are intended for the approximation of internal tissues and ligation of vessels. The PROXIMATE™ Reloadable Vascular Linear Stapler (product code TX30V) delivers three staggered rows of titanium staples and creates a 30 mm staple line.

The instruments may be reloaded during a single procedure. Do not reload the instruments more than seven times for a maximum of eight firings per instrument.

This device is intended to be used by healthcare professionals trained and experienced in surgical procedures involving the use of cutters for transection and resection.

A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>

One or more non-implantable components of this device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:		
	Material Present	Residual Risk
	Cobalt: CAS No. 7440-48-4 EC No. 231-158-0	<i>Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.</i>
The footnote adjacent to the above symbol on the packaging refers to the relevant Material Present		

MRI Safety Information

Non-clinical testing has demonstrated the staples made of titanium (Ti6Al4V) alloy are MR Conditional. A patient with the implanted staples can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 11,100 gauss/cm (111 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1.1 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the staples are expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the staple line when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

The product codes for the PROXIMATE™ Reloadable Vascular Linear Stapler and Reloads are as follows:

Instrument	Reload	Description	Number of Staples	Reload Color
TX30V	XR30V	Vascular	23	White

Note: The white (vascular) reload is dedicated to the 30 mm Vascular Linear Stapler and is not interchangeable with the blue and green reloads (sold separately).

Note: Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may require modification of the normally selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.

Illustration and Nomenclature (Illustration 1)

1. Anvil
2. Jaws
3. Release Button
4. Handle
5. Closing Trigger
6. Firing Trigger
7. Metal Housing
8. Staple Retaining Cap
9. Reload
10. Retaining Pin
11. Manual Pin Tab
12. Drivers

Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

1. Visually inspect all packaging for damage including a careful inspection of all sterile barrier systems for breaches in package integrity immediately prior to use.
Caution: Do not use if the product sterile barrier system or its packaging is compromised.
2. Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
3. Remove the staple retaining cap from the instrument. Discard the staple retaining cap (Illustration 2).
4. Position tissue to be stapled within the jaws of the instrument.

Note: Any tissue covering the hole in the anvil will be pierced by the retaining pin.

5. Squeeze the closing trigger until a click is heard. The instrument is in an intermediate position, the pin is fully seated in the anvil capturing the tissue (tissue not shown in illustration), and the jaws are partially open. Reposition tissue within the instrument if desired (Illustration 3).
Note: The retaining pin automatically seats into the anvil when operating the closing trigger. If desired, the retaining pin may be manually seated while the jaws remain fully open. Push the manual pin tab toward the anvil after positioning the tissue to be stapled within the jaws of the instrument.
Note: Inspect for the unwanted presence of clips, instruments, or other hard objects between the surfaces of the anvil and the staple cartridge (reload). Closing or firing over these items can damage the instrument and result in poor staple formation.
6. Squeeze the closing trigger and the handle fully together until a second click is heard. The closing trigger is now latched to the handle and the jaws are clamped onto the tissue (tissue not shown in illustration), which is ready to be stapled. The firing trigger will simultaneously move down to the ready-to-fire position (Illustration 4).
Note: Continue to grasp and manipulate the instrument using the closing trigger until ready to fire the instrument. Do not grasp the firing trigger before the instrument is to be fired.
Note: If the tissue needs to be repositioned within the instrument before stapling, open the jaws by pushing the release button and slowly releasing the grasp of the closing trigger. The closing trigger will return to the fully open position and the jaws will release the tissue. Tissue can now be repositioned.
Note: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple formation with the loss of staple line integrity and subsequent leakage, disruption, or poor healing. In addition, instrument damage or failure may result.
Caution: Unusually high closure force is a warning to open the instrument and inspect for tissue anomalies and hard objects or consider replacing the instrument.
7. Before firing, check to ensure the retaining pin is seated in the anvil. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.
Note: Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching, or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.
Note: When dividing vascular structures, it is mandatory to adhere to the basic surgical principle of proximal and distal control. Failure to fire the device, or incomplete firing of the device, and then proceeding to divide the vessel may lead to catastrophic bleeding. For this reason, inspect the staple line by releasing the device prior to dividing the vessel. An alternative approach is to place a vascular clamp across the vessel prior to dividing.
8. Fire the instrument by pulling the firing trigger back completely against the closing trigger until a click is heard, signifying the staples are fully formed (Illustration 5).
Note: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.
Note: The instruments and reloads are designed with a lock-out feature which prevents firing if a used reload or no reload is in the instrument. If the firing trigger does not pull back completely against the closing trigger, open the instrument as described in step 9. Replace the reload with a new reload.
Note: Make sure that the release button is not pressed during firing as proper formation of staples may be compromised.
9. Prior to opening the instrument and releasing the tissue, the edge of the reload or the side of the anvil may be used as a guide to transect tissue (not vessels) or to excise tissue (tissue not shown in illustration), which is protruding through the jaws. This aids in cutting at a proper distance from the staple line (Illustration 6).
Caution: For vascular applications, open the instrument and examine the integrity of the staple line before cutting.

10. Open the jaws by pushing the release button and releasing the grasp of the closing trigger (Illustration 7). The triggers and jaws will fully open, releasing the tissue. Remove the instrument.
Caution: Examine the staple line for hemostasis/pneumostasis and proper staple formation. If hemostasis/pneumostasis is not present, appropriate techniques should be used to achieve hemostasis/pneumostasis.

Reloading the Linear Stapler

1. Visually inspect all packaging for damage including a careful inspection of all sterile barrier systems for breaches in package integrity immediately prior to use.
Caution: Do not use if the product sterile barrier system or its packaging is compromised.
2. Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
3. Depress the release button to ensure the instrument is in the open position and the retaining pin is fully retracted into the reload.
4. Remove the used reload from the instrument. Grasp the top of the reload and lift upward, unsnapping the reload from the jaws. Properly discard the used reload (Illustration 8).
Caution: Inspect the instrument's anvil and jaw's surfaces after rinsing, in sterile solution, and wiping away any unused staples or debris.
5. Examine the new reload for the presence of a staple retaining cap. Remove the staple retaining cap by sliding it off the reload. If the retaining cap is not in place, discard the reload (Illustration 9).
6. Insert the new reload into the metal housing and snap into position. The tracks on each side of the reload should be used as guides to align the reload within the jaws of the instrument. When the reload is properly aligned, push the reload into the instrument until it is fully seated and a click is heard. Check that the reload is held firmly within the jaws (Illustration 10).
7. After reloading, observe the surface of the new reload. If colored drivers protrude out of the reload, replace with another reload. The linear stapler is now reloaded and ready for use.

Warnings and Precautions

- Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser and electro-surgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to other medical instruments.
- Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may require modification of the normally selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.
- Before firing, check to ensure the retaining pin is seated in the anvil. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.
- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis/pneumostasis, or disruption of staple line.
- Do not use if the product sterile barrier system or its packaging is compromised.
- Do not fire the instrument unless the closing trigger is properly latched against the handle.
- The firing trigger must be pulled back completely against the closing trigger to properly fire the instrument. Do not partially fire the instrument.

- The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.
- Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple formation with the loss of staple line integrity and subsequent leakage, disruption, or poor healing. In addition, instrument damage or failure may result.
- When dividing vascular structures, it is mandatory to adhere to the basic surgical principle of proximal and distal control. Failure to fire the device, or incomplete firing of the device, and then proceeding to divide the vessel may lead to catastrophic bleeding. For this reason, inspect the staple line by releasing the device prior to dividing the vessel. An alternative approach is to place a vascular clamp across the vessel prior to dividing.
- Examine the staple lines for hemostasis/pneumostasis and complete closure. If hemostasis/pneumostasis is not present, appropriate techniques should be used to achieve hemostasis/pneumostasis.
- The PROXIMATE™ Reloadable Vascular Linear Staplers may be reloaded during a single procedure. Do not reload the instruments more than seven times for a total of eight firings per instrument.
- Before closing the instrument on tissue, observe the surface of the new reload where the staples will emerge. If colored drivers protrude out of the reload, replace with another reload.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination. Disposal of all medical device instruments and accessories should be performed according to local requirements and regulations.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse, reprocessing, or resterilization of single-use devices may create a risk of contamination or cause infection or cross-infection, including but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the country in which the user and/or patient is established.

Storage and Transportation Conditions

Keep Dry
Keep Away from Heat

How Supplied

The PROXIMATE™ Reloadable Vascular Linear Staplers are supplied sterile and preloaded for single patient use. Discard after use.

The reloads for the PROXIMATE™ Reloadable Vascular Linear Staplers are supplied sterile for single patient use. Discard after use.

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©Ethicon Endo-Surgery, Inc. 2022
Johnson & Johnson Medical GmbH
Ethicon Endo-Surgery, LLC