	Registration document TECHNICAL SHEET	ENR 05/25 – Ind 3 – 11/06/24	
		FT-KIPIC	Rev: 4
		Date: 12/03/2025	P 1 / 4

DEVICE: KIPIC HYPODERMIC NEEDLES FOR SINGLE USE

1. ADMINISTRATIVE INFORMATION ON THE COMPANY		
1.1	Name: Aesthetic Group	
1.2	Address: Z.A. La Gobette 60540 Puisseux-Le-Hauberger France	Phone: +33 03 44 74 19 95 Fax: +33 03 44 74 18 94 E-mail: info@aestheticgroup.fr Website: www.aestheticgroup.fr
1.3	Vigilance and PRRC contact details: Valérie Boquet	Phone: +33 03 44 74 19 95 Fax: +33 03 44 74 18 94 E-mail: info@aestheticgroup.fr

2. INFORMATIONS ABOUT THE DEVICE OR EQUIPMENT

2.1 Name: Hypodermic needles for single use

2.2 Tradename: Kipic™

2.3 Nomenclature code:

UMNDS: 12-745

GMDN: 59230

HS Code: 90183210

2.4 LPPR code (if applicable): N/A

2.5 Class of the device:

Disposable Needles (Hypodermic Needles) is an invasive device intended for transient use. According to rule 6 of Annex VIII of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL, Disposable Needles (Hypodermic Needles) belongs to class II a medical device.

Number of the notified body: 0123

Date of first entry on the market in the EU: N/A

MD Manufacturer: Zhejiang kindly medical devices Co., Ltd.

Add: No.758, 5th Binhai Road, Binhai Industrial Park, Longwan District, 325025 Wenzhou, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA


SRN: CN-MF-000007594

BASIC UDI : 69230334202002a00401LY

2.6 Description of the device:

The Kipic™ disposable sterile needles are designed for hypodermic injections especially for infusion of drugs and blood collection, treat fat deposits and cellulite, alopecia, skin rejuvenation, pain management, osteoarthritis

They are available in different lengths and diameters. The needles consist of a straight cannula with bevel and a Luer Lock base.

	Registration document		ENR 05/25 – Ind 3 – 11/06/24	
	TECHNICAL SHEET		FT-KIPIC	Rev: 4
			Date: 12/03/2025	P 2 / 4

Procedure pack: No ☒ Yes

2.7 Catalogue references:

Item code	Color code	External diameter (mm)	Diameter (Gauge)	Length (mm)
KIPIC21G80	Green	0,82	21	80
KIPIC25G42	Orange	0,50	25	42
KIPIC27G04	Grey	0,40	27	4
KIPIC27G13	Grey	0,40	27	13
KIPIC27G42	Grey	0,40	27	42
KIPIC30G04	Yellow	0,30	30	4
KIPIC30G13	Yellow	0,30	30	13
KIPIC30G25	Yellow	0,30	30	25
KIPIC32G04	Green	0,23	32	4
KIPIC32G08	Green	0,23	32	8
KIPIC32G13	Green	0,23	32	13
KIPIC33G04	Black	0,20	33	4
KIPIC33G08	Black	0,20	33	8
KIPIC33G13	Black	0,20	33	13
KIPIC34G04	Orange	0,18	34	8

Packaging :

UCD (Control Unit): Blister of 1 piece

CDT (Multiple of the UCD): Boxes of 100 pieces

All the needles are equipped with a polypropylene transparent sterility protection tube. This tube protects the bevel of the needle so that it is not blunt and is suitable for ethylene oxide sterilization. The needles are packaged in an individual medical grade bag, which guarantees the protection of the needle and the maintenance of the sterile condition.

Labelling: other document

2.8 Composition of the device and accessories:


COMPONENTS	MATERIALS	TYPE OF CONTACT
Needle tube	Medical grade 304 stainless steel	Direct
Protective cap	Polypropylene	Indirect
Needle hub	Polypropylene	Direct

Latex:	Presence	Absence
Phthalates (DHP):	Presence	Absence
Product of animal or biological origin:	Presence	Absence
Bisphenol A	Presence	Absence

Associated devices and accessories: According to needle variants, Disposable Needles (hypodermic needles) should be used with syringe. Requirements for Disposable Needles (hypodermic needles) have been specified in ISO 7864.

3. STERILIZATION PROCESS

MD sterile: YES ☒ NO

	Registration document TECHNICAL SHEET	ENR 05/25 – Ind 3 – 11/06/24	
		FT-KIPIC	Rev: 4
		Date: 12/03/2025	P 3 / 4

Sterilization method of the device: Ethylene Oxide

Shelf life: 5 years

4. PRESERVATION AND STORAGE

Store in a dry environment and keep away from sunlight.

5. SAFETY OF USE

5.1 Technical safety:

The MRI compatibility of our needles has not been tested

Needle is X-ray detectable

Like all stainless steel needles, steel could be altered by hydrochloric acid

5.2 Biological safety:

Be careful to follow the usual aseptic procedures

After use, handle and dispose of with other contaminated devices

6. INSTRUCTIONS FOR USE

6.1 Instructions for use:

1. Select appropriate matching device with standard luer and gauge as per actual requirement.
2. Check package integrity before use.
3. Check the integrity of the package before use; if the primary package is opened or damaged or in case there is foreign matters such as hair, insects, do not use the product.
4. Open the package and take out the product from the package. Take care to prevent the product falling to the ground.
5. Inspect product, do not use the product in case of abnormal product conditions including protective cap detachment, inverted or unbonded needle tube.
6. Assembled it with the standard luer of disposable syringe.
7. Perform hypodermic injection on the human body according to clinical nursing science operating method.
8. Dispose the product after use in conformity with the institution's regulatory for contaminated sharps and according to valid national regulations for such products. Take care to avoid accidental needle stick injury during sharp disposal.

6.2 Intended purpose:

Disposable Needles (Hypodermic Needles) is intended to inject or withdraw fluids from primarily the human body. The product is intended to be used with disposable syringe, the needle shall be used immediately after medication filling.

6.3 Intended users:

Healthcare professionals


6.4 Intended patient population:

Not limited, mainly pediatrics and adults

6.5 Indications:

Patients who need hypodermic injection therapy

6.6 Using precaution and warnings:

	Registration document		ENR 05/25 – Ind 3 – 11/06/24	
	TECHNICAL SHEET		FT-KIPIC	Rev: 4
			Date: 12/03/2025	P 4 / 4

- The product must not be reused.
- Do not resterilize
- Entrust the operation and use of the product and accessories only to persons with the required training, knowledge or experience.
- Physicians administering drug therapy should be aware of the chemicals / medications injected, dosages, particular side effects and potential interactions.
- Use the product only for the intended purpose, see Scope of application / see Contraindications.
- Before each use, make a visual examination of the product: absence of loose, bent, broken, cracked, worn and broken parts.
- Never use a damaged or defective product. Dispose of any damaged product immediately.
- Check the integrity of the primary packaging before use. Never use this device when the primary packaging is damaged.
- Do not use beyond the expiry date
- Be very careful to apply the usual aseptic procedures due to the risk of infection
- Do not exert excessive force that may cause the needle to break
- Solutions based on hydrochloric acid can alter the steel
- Check the integrity of the needle during removal
- For single use only. Discard after use
- After use, handle and dispose of contaminated devices - Do not recap
- Use immediately after opening the packaging and avoid polluting

6.7 Contra-indications:

N/A

7. ADDITIONAL PRODUCT INFORMATION

Biocompatibility

Comply with the requirements of the standard for biocompatibility of medical devices.

8. ANNEXES

- EC Certificate
- EU Declaration of conformity
- IFU
- Warnings

9. PICTURES

