Nasopharyngeal Airway with Adjustable Flange Instructions for use	en
Nasopharyngealtubus mit verstellbarem Haltering Gebrauchsanweisung	de
Sonde naso-pharyngée avec collerette réglable Consignes d'utilisation	fr
Via aérea nasofaringea con pestaña adaptable Instrucciones de uso	es
Tubo nasofaringeo con flangia regolabile Istruzioni per l'uso	it
Tubo nasofaríngeo com flange ajustável Instruções de utilização	pt
Nasofaryngeal luftvägstub med justerbar fläns Bruksanvisning	sv
Nasofaryngeale tube met instelbare flens Gebruiksaanwijzing	nl
Nasofaryngeal tube med justerbar flange Brugervejledning	da
Ayarlanabilir flanşlı Nazofaringeal Havayolu Kullanım talimatı	tr
Rurka nosowo-gardłowa z regulowanym kolnierzem Instrukcja obsługi	pl
可調フランジ付き鼻咽頭エアウェイ 取扱説明書	ja
帶可調法兰的鼻咽通气导管 使用说明	zh
조절가능 플랜지형 비인두용 에어웨이 사용 설명서	ko

3	Instructions for use	en	
5	Gebrauchsanweisung	de	
7	Consignes d'utilisation	fr	
9	Instrucciones de uso	es	
11	Istruzioni per l'uso	it	
13	Instruções de utilização	pt	
15	Bruksanvisning	sv	
17	Gebruiksaanwijzing	nl	
19	Brugervejledning	da	
21	Kullanım talimatı	tr	
23	Instrukcja obsługi	pl	
25	取扱説明書	ja	
27	使用说明	zh	
29	사용 설명서	ko	

PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE USE!

Product Description

Nasopharyngeal Airway with Adjustable Flange consists of an adjustable flange and tube.

Nasopharyngeal Airway constructed of green, Mediprene material with product number beginning with 185420 is **not** made with natural rubber latex.

Nasopharyngeal Airway constructed of yellow, Latex material with product number beginning with 185200 contains natural rubber latex which may cause allergic reactions and consequently should not be used in patients with natural rubber latex allergy!

Intended Use

To provide a patent airway for a patient with a blockage or potential blockage of the nasal cavity or nasopharynx.

Contraindications

Latex allergy (for product contains natural rubber latex only), traumatic brain injury or central facial fractures, basilar skull or cribriform fractures, anticoagulant therapy, coagulation disorders, transphenoidal surgery, and new postoperative rhinoplasty or septoplasty because it can cause tissue trauma or damage the newly altered structural integrity of the surgical site.

Direction for Use

1) Lubricate the airway with a small amount of water-based lubricant avoiding the part of tube where the adjustable flange is located.

Caution: Excessive use of lubricant can cause the adjustable flange to slip along the tube easily and may dislocate the position of the tube in situ.

2) Gently insert the airway into the nostril. Do not force. Gently rotate the airway to ease entry. If insertion is not possible, apply the tube via the contralateral nostril.

3) Adjust the adjustable flange after proper positioning of the tube. Proper positioning is achieved when the tip of the tube is visible at the level of the uvula.

Caution: Do not use excessive force as it may separate the adjustable flange from the cannula.

4) If a suction catheter is used through the central lumen of the nasopharyngeal airway, ensure correct size and lubricate the outer surface of the suction catheter prior to inserting it into the airway for mucous clearance/removing body fluids purpose. Grasp the tube's funnel with one hand while inserting the suction tube with the other hand.

Caution: Avoid using strong force when introducing the suction catheter into the airway.

5) If face mask is used in conjunction with the nasopharyngeal airway, ensure fitting of face mask does not press against or push the airway into the nostril.

 Warning

Some products contain natural rubber latex which may cause allergic reactions. Please refer product label.

CAUTION:

- Do not force the system through the nose, should significant resistance exist.
- Do not attempt insertion where lesions or infection in the nose is observed.

Re-processing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected, or sterilized.

Storage Instruction

Keep away from sunlight and keep dry. Do not use if the product sterilisation barrier or its packaging is compromised.

**Caution: Federal Law restricts this device to sale by or on the order of a physician.
(For USA and Canada only)**

Explanation of important symbols and markings on the product label.

SIZE	Product size
Ch.	1 Ch. = 1/3mm
Fr.	1 Fr. = 1/3mm
O.D.	Outer diameter
Qty	Quantity, Number of items
	Not made with natural rubber latex
	This product contains natural rubber latex
	Do not use if package is damaged
	"Green Dot" Package recycling

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