Safe Ear Currettes™

INSTRUCTIONS FOR USE	
MD	CE
Ⓡ	EC REP Obelis sa Boulevard Général Wahis 53 1030 Brussels, BELGIUM Tel: +(32) 2.732.59.54 Fax: +(32) 2.732.60.03 E-Mail: mail@obelis.net
Rx ONLY	Ⓡ Bionix LLC 1670 Indian Wood Circle Maumee, Ohio 43537, United States Phone: +1 419.727.8421 Fax: +1 419.727.8426 www.Bionix.com
MADE IN THE U.S.A.	

Importer, Importeur, Importör, Importador, Importateur, Importeur, Importator, Dovezce, Maahantuaja, Importatore, Έισαγωγέας, Importuotojas:

MedEnvoy
Primes Magrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

This manual applies to the **REF** numbers indicated below:

#1222 - Green MicroLoop®
#3444 - Red AngleLoop®
#4888 - Blue InfantScoop®
#9555 - White FlexLoop®

#2999 - Orange ControLoop®
#4111 - Purple VersaLoop®
#6333 - Yellow CeraSpoon®
#5777 - Variety Pack

	©2022 Bionix LLC 1670 Indian Wood Circle Maumee, Ohio 43537 United States Phone: 800.551.7096 Fax: 800.455.5678 Email: bionix@bionix.com www.bionix.com RM95-0005 2022/05/11 Rev. F
bionix ™	
ENGLISH	
INSTRUCTIONS FOR CERUMEN REMOVAL	

Intended use: The intended use of the Safe Ear Curette is to remove cerumen and/or foreign bodies from a patient’s ear canal without causing injury.

Indications for use: The Safe Ear Curette is indicated for use in cases of impacted cerumen and/or foreign bodies occluding the ear canal.

Position the patient comfortably on the examination table. Smaller children may be restrained on their parent’s lap or shoulder during the procedure.

Using an otoscope, examine the ear canal and tympanic membrane, noting the presence and location of any cerumen in the ear canal.

Select an ear curette of the appropriate size for the patient:

White FlexLoop® for children and adults, general purpose; **Blue InfantScoop®** for smaller infants and children; **Green MicroLoop®** for the harder cerumen of older children and adults; **Yellow CeraSpoon®** for fast removal of extensive cerumen **Red AngleLoop®** for hard to reach cerumen; **Purple VersaLoop®** for gentle precision; and the **Orange ControLoop®** for precise control and direction during curettage.

Straighten the ear canal by gently pulling the ear up and away from the patient. Gently use the curette to loosen and draw cerumen from the ear canal. Use an otoscope to frequently re-examine the ear canal during the procedure as cerumen may shift position or move during curettage.

Note the presence and location of the remaining cerumen, if any. Repeat the curettage procedure until sufficient cerumen is removed to allow adequate examination of the tympanic membrane.

CAUTION: Discontinue curettage immediately if bleeding, irritation, or other trauma to the ear canal or tympanic membrane occurs.

WARNING: Cross-contamination risk. Do not reuse disposable curette tips as this may spread contamination from one patient to another patient.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Clinical Benefits: Indicated in Intended Use.

Contraindications: There are no absolute contraindications for cerumen removal.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

MD	Medical Device	Ⓡ	Do Not Use if Package is Damaged
Ⓡ	Manufacturer	Ⓡ	Do Not Reuse Single-Use Device
Ⓢ	Date of manufacture	Ⓢ	Single patient-multiple use
Ⓜ	Use-By Date	Ⓜ	Consult Instructions For Use (IFU)
LOT	Lot Code	!	Caution
REF	Reorder Number	!	Warning
SN	Serial number	EC REP	Authorized representative in the European Community
STERILE R	Sterilized using irradiation	CE	European Conformity
Ⓡ	Do Not Resterilize	Rx ONLY	Prescription Only or "For Use by or on the order of a licensed medical professional"
⚠	Non-sterile	Ⓢ	Temperature limitation

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