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| | |
|--|----|
| Product Description and Intended Use | 30 |
| Warnings & Cautions, Contraindications | 31 |
| Principles of Operation | 32 |
| Assembly View | 35 |
| Replacement Parts | 36 |
| Operating Instruction | 38 |
| Cleaning, Disinfection and Sterilization | 40 |
| Storage | 47 |
| Testing the Manual Resuscitator | 48 |
| Specifications & Performance Characteristics | 50 |
| The PEEP Valve (Accessory Series) | 53 |

EN Product Description and Intended Use

The manual resuscitator is designed for use as an aid to artificial respiration and cardiopulmonary resuscitation. The manual resuscitator is used to ventilate apneic patients and to support breathing and/or supply oxygen to spontaneously breathing patients. The resuscitator has different designs for adults, children and infants and is adapted to the different compressing frequencies. Different sizes are available according to the patient's oxygen requirements.

An optional pressure relief valve can be provided. When the pressure inside the bag exceeds 60 cm H₂O for adults and 40 cm H₂O for children, the pressure relief valve will automatically vent the supplied ventilation to the atmosphere to protect the lungs from a barotrauma.

This product is intended for use by qualified medical or emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques. This manual provides all section assembly drawings, all parts numbers and information on cleaning procedures. Before using this manual resuscitator, please read this manual carefully to accomplish the best effect, and pay attention to all safety warnings.

Silicone Manual Resuscitator

This product is a medical device.

It is made of stable silicone, which is extremely temperature-resistant (up to a maximum of 134 °C).

EN Warnings & Cautions, Contraindications

WARNINGS

1. Do not use the manual resuscitator in toxic atmospheres.
2. If no additional oxygen is being administered, remove the oxygen reservoir. Failure to do so will affect the refill rate and maximum frequency capabilities.
3. Do not administer additional oxygen in the presence of open flames.
4. Do not use any oil/grease or other hydrocarbon-based substances on any part of the manual resuscitator. Additional oxygen administered under pressure can combine with hydrocarbons and cause explosions.
5. This device is intended for use by qualified medical and emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques.
6. Proficiency in the assembly, disassembly and use of this product should be demonstrated before use on a patient.
7. Always check the product in accordance with this manual after cleaning and sterilization or replacement of parts.
8. Monitor airway pressure with a manometer when ventilating a patient.
9. Only qualified personnel trained in the use of positive end expiratory pressure (PEEP) should administer PEEP with this product.
10. Verify the PEEP level and the function of the resuscitator before use on the patient.

CAUTIONS

1. If deactivating the pressure relief valve, it must be ensured that the pressure in the patient's airways does NOT become too high.
2. Do not attempt to disassemble the pressure relief valve. Disassembly will damage the components.
3. Before use, clean and sterilize the resuscitator according to your institution's validated procedure. Then check the product functions correctly as specified in this manual.


CONTRAINDICATIONS

Not suitable for patients who have allergies to the material.

EN Principles of Operation

Assembly View (I)

The resuscitator consists of 4 components

Assembly View (I)


- | | |
|---|---|
| (A) Patient valve (consists of E, F, H) | (E) silicone membrane (duckbill design) |
| (B) Silicone bag | (F) expiratory/inspiratory port |
| (C) Intake valve | (G) gas inlet |
| (D) Oxygen reservoir | (H) PEEP valve adapter |

If no additional oxygen is being administered, the oxygen reservoir should be removed.

Principle Diagram (II) - Inspiration

When you press the resuscitator, positive pressure is generated and the intake valve (C) closes. The air inside the bag presses on the duckbill valve (E), thereby blocking the expiratory/inspiratory port (F) and peep valve adapter (H); air is therefore transported to the patient's lungs via the silicone membrane.

If oxygen has to be supplied, it should be connected to gas inlet (G). In this way, the oxygen fills the reservoir through the intake valve, passes into the silicone bag by means of the inspiratory action and then flows into the patient when the bag is squeezed.


Principle Diagram (II) - Inspiration

EN Principles of Operation

Principle Diagram (III) - Exhalation

When the silicone bag (B) is released, the silicone membrane remains closed, with the result that the exhaled air flows out via the expiratory/inspiratory port (F). At the same time, the intake valve (C) is opened by the expiratory pressure created when the bag is released and air flows into the bag through the top of the intake valve. At the same time, oxygen flows from the oxygen reservoir into the bag until the bag has resumed its original shape that it had before compression.

Principle Diagram (III) - Exhalation


To prevent an excessive oxygen flow rate and a low compression frequency causing excessively high pressure inside the bag and reservoir, the intake valve (C) is specially designed to release the excessive air to maintain a low oxygen supply rate and ensure the patient's safety.

EN Principles of Operation

Principle Diagram (IV)


The resuscitators for infants and children are equipped with pressure relief valves and automatically adjust the pressure in the lungs to $40 \text{ cm H}_2\text{O} \pm 5 \text{ cm H}_2\text{O}$. If the pressure is higher than this, the pressure relief valve opens for safety reasons and releases the pressure.

Pressure relief valve (pop-off valve)


Activating the pressure relief valve

When higher inspiratory pressure is required, the pressure relief valve can be deactivated by pressing down and turning the valve as shown.

Deactivating the Pressure Relief Valve

EN Assembly View

EN Replacement Parts

| REF | Resuscitator set | Area of application | Pressure relief valve | Volume | Material |
|---------------|---|---------------------|------------------------|---------|------------------------|
| 51.500.00.458 | Resuscitator with pop-off-valve and PEEP valve adaptor | Infant | 40 cm H ₂ O | 280 mL | Silicone, PSU |
| 51.500.00.454 | Resuscitator with pop-off-valve and PEEP valve adaptor | Child | 40 cm H ₂ O | 500 mL | Silicone, PSU |
| 51.500.00.453 | Resuscitator with pop-off valve and PEEP valve adaptor | Adult | 60 cm H ₂ O | 1600 mL | Silicone, PSU |
| 51.500.00.452 | Resuscitator with PEEP valve adaptor | Adult | – | 1600 mL | Silicone, PSU |
| 51.500.00.300 | Resuscitator with PEEP valve adaptor, two masks size 4 + 5 | Adult | – | 1600 mL | Silicone, PSU, PC |
| 51.500.00.301 | Resuscitator with PEEP valve adaptor, two masks size 4 + 5, Oxygen reservoir: 2500 mL | Adult | – | 1600 mL | Silicone, PSU, PC, PVC |
| 51.500.00.302 | Resuscitator with PEEP valve adaptor, two masks size 2 + 3, Oxygen reservoir: 2500 mL | Child | 40 cm H ₂ O | 500 mL | Silicone, PSU, PC, PVC |
| 51.500.00.303 | Resuscitator with PEEP valve adaptor, two masks size 0 + 1, Oxygen reservoir: 600 mL | Infant | 40 cm H ₂ O | 280 mL | Silicone, PSU, PC, PVC |
| 51.500.00.100 | Resuscitator with PEEP valve adaptor, two masks size 0 + 1 | Infant | 40 cm H ₂ O | 280 mL | Silicone, PSU, PC |
| 51.500.00.200 | Resuscitator with PEEP valve adaptor, two masks size 2 + 3 | Child | 40 cm H ₂ O | 500 mL | Silicone, PSU, PC |

PSU – Polysulphone
 PC – Polycarbonate
 PVC – Polyvinyl chloride

EN Replacement Parts

| REF | Accessories | Area of application | | Volume | Material |
|---------------|-------------------------------|---------------------|--|---------|----------|
| 51.500.00.538 | Oxygen reservoir (single use) | Adult Child | | 2500 mL | PC, PVC |
| 51.500.00.540 | Oxygen reservoir (single use) | Infant | | 600 mL | PC, PVC |

| REF | Accessories | Ø OD/ID | Max. PEEP level | Min. PEEP level | Material |
|---------------|-------------------------|-------------|--------------------------|-------------------------|-------------------------------|
| 51.500.00.043 | PEEP valve (single use) | 22 mm/15 mm | 10.0 cm H ₂ O | 2.0 cm H ₂ O | Stainless steel, PC, Silicone |
| 51.500.00.045 | PEEP valve (single use) | 22 mm/15 mm | 20.0 cm H ₂ O | 5.0 cm H ₂ O | Stainless steel, PC, Silicone |

| REF | Accessories | Area of application | | Material |
|---------------|---|----------------------|--|---------------|
| 51.500.00.528 | Intake valve | Adult, Child, Infant | | PSU, Silicone |
| 51.500.00.404 | Pop-off patient valve, PEEP valve adaptor | Adult | | PSU, Silicone |
| 51.500.00.400 | Patient valve with PEEP-valve adaptor | Adult | | PSU, Silicone |
| 51.500.00.409 | Silicone membrane for patient valves, duckbill design | Adult | | Silicone |

| REF | Ventilation masks | Area of application | | Material |
|---------------|-----------------------|---------------------|--|--------------|
| 51.500.00.000 | Silicone mask, size 0 | Infant | | Silicone |
| 51.500.00.001 | Silicone mask, size 1 | Infant | | Silicone |
| 51.500.00.002 | Silicone mask, size 2 | Child | | Silicone |
| 51.500.00.003 | Silicone mask, size 3 | Child | | Silicone, PC |
| 51.500.00.004 | Silicone mask, size 4 | Adult | | Silicone, PC |
| 51.500.00.005 | Silicone mask, size 5 | Adult | | Silicone, PC |

PSU – Polysulphone
 PC – Polycarbonate
 PVC – Polyvinyl chloride

1. The patient is in the supine position – pull his/her chin up to keep the airway and the oral cavity in alignment. This will make it easier to breathe.

2. Remove all material from the mouth and throat.

3. Insert the oropharyngeal tube into the mouth to keep it open and to prevent the tongue from blocking the airway. (A mouth gag can be used to open the mouth). Select the oropharyngeal tube according to the size of the oral cavity.

4. The emergency personnel should stand behind the patient's head, extend the head back and pull the chin up and toward the emergency personnel.


Note

If the patient already has an airway tube inserted or has undergone a tracheostomy, no need to use the mask but connect the patient valve connector to the tube, then proceed as usual.

5. Cover the patient's mouth and nose with the mask and press on the mask to maintain an airtight seal with the patient's face.

6. Use the other hand to squeeze the resuscitator with a sufficient inhalation and exhalation frequency. (Adult: 12-16 times, child: 14-20 times, infant: 35-40 times)

7. The emergency personnel must check that the patient is being correctly ventilated.

- Check the patient's breathing movements.
- Check the color of the patient's lips and face through the transparent part of the mask.
- Check that the patient valve is functioning correctly through the transparent housing.
- Check whether the inside of the mask is foggy during exhalation.

Information regarding Cleaning, Disinfection and Sterilization

- The resuscitator and accessories must be disassembled as shown on page 35. Do not disassemble the parts further than shown.
- For parts exposed to expiratory gases, cleaning, disinfection and sterilization must be performed after each patient. For parts not exposed to patient expiratory gases, cleaning, disinfection and sterilization must be performed regularly, as needed, to remove dust etc. If the resuscitator has been used for patients/in environments with infectious diseases, cleaning, disinfection and sterilization must be carried out for the entire resuscitator set.

Risk Assessment and Classification of Medical Devices prior to Processing

The nature and scope of processing are determined according to the use of the medical device. For this reason, the operator is responsible for correctly classifying the medical devices and, therefore, for defining the nature and scope of processing. On the basis of this user-dependent classification, the operator is able to define which of the processing methods listed in this processing instruction must be performed.

Note the manufacturer's information about material compatibility for cleaning, disinfection and sterilization.

INSTRUCTIONS

Storage and Transport

- It is recommended that the contaminated medical devices be transported in a sealed container.
- It is recommended that the medical devices be processed as soon as possible after use and no more than 2 hours after use. Intermediate storage of used medical devices with contamination in the form of blood traces, for example, may result in damage to the device.

Preparation

- Wear personal protective equipment (sturdy gloves, water-repellent protective gown, face mask or safety goggles; also note the instructions of the disinfectant and cleaning agent manufacturer).
- Prior to cleaning and disinfection, some medical devices must be opened as much as possible and/or disassembled in accordance with the manufacturer's instructions for use.
- Do not remove the springs of the **pop-off valve!**

Pretreatment

Equipment for patient valve with integrated PEEP adapter and pop-off valve, intake valve, silicone membrane and silicone bag of resuscitator:

Washing gun (3.8 bar), tap water ($20 \pm 2^\circ\text{C}$) (at least drinking water grade)

1. Spray each medical device for at least 1 minute with a high-pressure washing gun in order to remove coarse surface contamination.

Equipment for silicone face masks:

Plastic brush (brush head length: 25 mm and 35 mm, bristle length: 5 mm and 10 mm), running water (at least drinking water grade), tub.

1. Fill tub with water. Place all masks in the tub with the water. All devices should be placed in such a way that they are completely under the surface of the water.
2. Remove coarse surface contamination on the medical devices by brushing them under the surface of the water with a plastic brush (an exposure time of 2 minutes under the water surface has been validated).

Cleaning: manual

Note: Coarse surface contamination on the medical devices must be removed prior to manual processing (see pretreatment).

Equipment for patient valve with integrated PEEP adapter and pop-off valve:

Multi-stage enzymatic cleaner (Sekusept MultiEnzyme P has been validated), **cleaning brushes:** brush head length: 70 mm, diameter: 15 mm, brush head length: 30 mm, diameter: 2.5 mm, brush head length: 13 mm, diameter: 1.6 mm, tap water/running water (at least drinking water grade) ($20 \pm 2^\circ\text{C}$)

Equipment for intake valve and silicone membrane:

Multi-stage enzymatic cleaner (Sekusept MultiEnzyme P has been validated), **cleaning brushes:** brush head length: 70 mm, diameter: 15 mm, brush head length: 30 mm, diameter: 2.5 mm, tap water/running water (at least drinking water grade) ($20 \pm 2^\circ\text{C}$).

Equipment for the silicone bag of the resuscitator:

Multi-stage enzymatic cleaner (Sekusept MultiEnzyme P has been validated), **cleaning brushes:** brush head length: 70 mm, diameter: 15 mm, tap water/running water (at least drinking water grade) ($20 \pm 2^\circ\text{C}$).

Equipment for silicone face masks:

Multi-stage enzymatic cleaner (Sekusept MultiEnzyme P has been validated),
plastic brushes: brush head length: 25 mm and 35 mm, bristle length 5 mm and 10 mm, tap water/running water (at least drinking water grade) (20 ± 2°C)

1. Prepare the cleaning solution according to the manufacturer's instructions (Sekusept MultiEnzyme P 2% has been validated).
2. Disassemble the medical devices and completely submerge them under the liquid level of the cleaning solution without creating air bubbles.
3. During the exposure time, brush the hard-to-reach areas of the submerged medical device with a soft brush (1 minute has been validated). Pay special attention to areas that are critically difficult to reach.
4. Exposure time in the cleaning solution at least 5 minutes (5 minutes has been validated).
5. Remove the medical devices from the cleaning solution and rinse each device thoroughly under running water (30 seconds has been validated).
6. Check for cleanliness; if contamination is still visible, repeat the previous steps.

Disinfection: manual with all components

Equipment for the resuscitator with all components:

Virucidal disinfectants (listed by VAH [German Association for Applied Hygiene]) for medical devices, e.g. based on peracetic acid (Ecolab Sekusept aktiv #3050570 has been validated), if at all possible, demineralized water (20 ± 2°C) free from facultative pathogenic microorganisms, tap water (20 ± 2°C), tub for disinfectant, lint-free disposable wipe, medical-grade compressed air.

For the following components brushes are needed.

Cleaning brushes for patient valve with integrated PEEP valve adapter and pop-off-valve (the following Interlock brushes have been validated:

brush head length: 70 mm, diameter: 15 mm,
 brush head length: 30 mm, diameter: 2.5 mm,
 brush head length: 13 mm, diameter: 1.6 mm).

Plastic brush for silicone face mask, brush head length: 25 mm and 35 mm, bristle length: 5 mm and 10 mm)

1. Prepare the disinfectant solution according to the manufacturer's instructions (Sekusept aktiv 2% has been validated).
2. Completely submerge the medical devices in the disinfectant solution and make sure that all accessible surfaces are wetted.
3. Applicable for patient valve with integrated PEEP adapter, pop-off valve as well as the silicone face masks: During the exposure time, brush the hard-to-reach areas of the submerged medical device with a soft brush for at least 1 minute (1 minute has been validated). Pay special attention to the areas that are critically difficult to reach.
4. Exposure time in the disinfectant solution according to the disinfectant manufacturer's instructions (Sekusept aktiv 2% for 5 minutes has been validated).
5. Remove the medical devices from the disinfectant solution and allow to drip dry.
6. Rinse the medical devices for 1 minute with demineralized water (1 minute has been validated).
7. Wipe with a lint-free disposable wipe or dry with medical-grade compressed air.

Cleaning and disinfection by machine

Note: Coarse surface contamination on the medical devices must be removed prior to machine processing (see pretreatment).

Equipment:

Washer-disinfector according to DIN EN ISO 15883-1+2 with thermal program (temperature 90°C to 95°C), cleaning agent: mild alkaline cleaning agent (Sekumatic MultiClean has been validated), medical-grade compressed air.

1. Place all individual parts of the medical devices in suitable strainer bowls and small-parts baskets or place them on the shelf rack in such a way that all inner and outer surfaces of the medical devices are cleaned and disinfected.
 - 1.1 Place the silicone membrane of the resuscitator in a small-parts strainer.
 - 1.2 For complete rinsing through the inner lumina, the resuscitators must be placed upright in a basket or fitted on the shelf rack. The fit must be checked for firm seating prior to starting the machine and at the end of the process (upright positioning in a basket has been validated).

EN Cleaning, Disinfection and Sterilization

2. Close the washer-disinfector and start the program; see the table below for the program sequence.

| Program steps | Water | Dosage | Time | Temperature |
|--------------------------|---------------------|--|---------|--|
| Prerinseing | Cold water | | 5 min. | |
| Dosing of cleaning agent | | According to the manufacturer's instructions | | According to the manufacturer's instructions |
| Cleaning | Demineralized water | | 10 min. | 55°C |
| Rinsing | Demineralized water | | 2 min. | |
| Disinfection | Demineralized water | | 3 min. | A ₀ value > 3000* (e.g.: 90°C, 5 min.) |
| Drying | | | 15 min. | Up to 120°C |

* Within their jurisdiction, some authorities may issue different implementing regulations (parameters for disinfection performance).

3. After the end of the program, remove all individual parts of the medical devices.
4. Check for dryness of the load and, if necessary, dry with medical-grade compressed air.
5. After removal from the washer-disinfector, perform a visual cleanliness inspection. If contamination is still visible, clean the medical devices again manually. The cleaned medical devices must then be processed by machine again.

EN Cleaning, Disinfection and Sterilization

Maintenance, Inspection and Testing

A visual inspection of all medical devices for cleanliness, integrity and functionality must be performed, if necessary using an illuminating magnifier (3-6 diopter).

All medical devices (including parts of medical devices that can be disassembled) must be checked for damage and wear.

Damaged medical devices may no longer be used and must be rejected.

Packaging

Equipment: Foil-paper packaging (steriCLIN item no. 3FKFS230110 and 3FKFS230112 have been validated), seam sealer (HAWO Type 880 DC-V has been validated).

1. Package silicone face mask, in disassembled condition.
2. Reinstall the silicone membrane into the patient valve.
3. A suitable method (sterile barrier system) must be used to package the medical devices. Packaging in accordance with DIN EN ISO 11607 (individual) or DIN 58953-9 (sets).

Individual: A sterile barrier system (for example, foil-paper packaging) according to DIN EN ISO 11607, which has been designated by the manufacturer as suitable for steam sterilization, must be used. The medical devices may be singly or doubly packaged. The packages must be large enough that the seal seam is not under tension (double packaging in foil-paper packaging has been validated).

Note: After the heat-sealing process, the seal seam must be visually checked for possible defects. If defects are found, the package must be opened and the medical device repackaged and resealed.

Sterilization

Machine: Sterilizer according to DIN EN 285 or small steam sterilizer according to DIN EN 13060, Class B method.

Method: Steam sterilization with fractionated pre-vacuum, 134°C, holding time at least 3 minutes or 132°C for at least 3 minutes (longer holding times are not possible).

EN Cleaning, Disinfection and Sterilization

1. Place the packaged devices into the sterilizer chamber.
2. Start the program.
3. After the end of the program, remove the devices and allow them to cool.
4. Then check the packages for any damage and moisture penetration. Rejected packages must be considered non-sterile. The medical devices must be repackaged and resterilized.

Storage

Storage (protected against dust and humidity as a minimum) and storage time has to be according to the specification and internal requirements of the hospital.

Additional Information

Only validated processes may be used for the processing of medical devices.

Visual Inspection through Package Assembly

1. After the cleaning, inspection and sterilization process, let the resuscitator parts dry at room temperature. Wait until all parts are dry.
2. After the cleaning, disinfection and sterilization process, carefully inspect all parts for damage or excessive wear and replace them if necessary. Some methods may cause discoloration of rubber parts but will not affect their lifetime. In the event of material wear, e.g. cracks, the affected parts should be replaced. Please contact your distributor for parts replacement.
3. Assemble the parts as per the assembly view.
4. After assembling the resuscitator, a label should be affixed to indicate the handling date.

Note: This product has been validated for 40 reprocessing cycles.

EN Storage

- For compact storage (e.g. in an emergency case), the inlet end at the adult type can be pushed halfway into the bag.
- Never store the resuscitator in a compressed or folded state.
- Never excessively squeeze the bag during storage. When the resuscitator is ready for use, it should not be stored in direct sunlight or in a heated environment.
- Storage temperature: +5°C to +40°C
- For long-term storage or transportation, the resuscitator should be stored in a sealed package. Do NOT store in the sun or a heated environment.

The resuscitator should be tested as follows:

1. Test before use
2. After cleaning and sterilization
3. After fitting new parts
4. Monthly check if the resuscitator is not used frequently

Equipment required:

Test lung, 0 - 60 cm H₂O manometer, flow meter, regulated gas supply, gas supply tubing.

Testing the Silicone Bag

1. Remove the patient valve and the oxygen reservoir.
2. First squeeze the silicone bag, afterwards manually close the opening, where you attached the Patient valve.
3. The bag should expand immediately and refill. If not, check whether the intake valve at the base of the silicone bag is correctly assembled.
4. First close the opening, where you attached the Patient valve, afterwards squeeze the silicone bag. The bag should not be easy to squeeze. However, if it is easy to squeeze, check whether the opening is manually correctly closed and the intake valve at the base of the silicone bag is correctly assembled.

Testing the Patient Valve

1. Connect the patient valve to the silicone bag. Connect the test lung to the outlet of the patient valve.
2. Squeeze and hold the bag. The silicone membrane (duckbill design) in the patient valve should open and the test lung should fill up.
If this is not the case, check the connection between the resuscitator and the test lung, and check that the patient valve is correctly assembled.
3. The silicone membrane (duckbill design) should close, and as the test lung deflates, gas should flow through the expiratory/inspiratory port.
If not, check whether the patient valve with the silicone membrane is correctly assembled.
4. Ventilate the test lung for a minimum of ten cycles to ensure that the resuscitator is functioning correctly. Inspiration should occur when the silicone bag is squeezed; expiration should occur when the bag is released. If not, check whether the patient valve with the silicone membrane is correctly assembled.

Functional Check of the Pressure Relief Valve (Infant and Child Resuscitator)

Connect a manometer (0-100 cm H₂O) to the patient outlet of the patient valve. Squeeze the bag. When the pressure relief valve is activated, the value measured by the manometer should be 35 - 45 cm H₂O. If not, check whether the patient valve is correctly assembled and that there are no leaks. If the pressure relief valve does not pass the test, it must be replaced. Do not attempt to repair the pressure relief valve.

Overall Resuscitator Function

1. Fully assemble the resuscitator (patient valve, silicone bag, intake valve and oxygen reservoir). Connect the resuscitator to an external gas source and connect a test lung to the patient valve with silicone membrane.
2. Set the external gas flow to 15 liters per minute for adults and children and to 10 liters per minute for newborns and infants.
3. Perform several ventilation cycles. The test lung should inflate during inspiration and deflate during expiration. Check for leaks at all joints and connections. Make sure the resuscitator fills properly and all valves are operating correctly.

EN Specifications & Performance Characteristics

Operating temperature: -18°C to 50°C

Connections

| | |
|----------------------|----------------------|
| Patient valve | ID: 15 mm; OD: 22 mm |
| Silicone bag inlet | ID: 23 mm |
| Intake valve port | OD: 25 mm |
| Additional gas inlet | OD: 6 mm |

Dead space

| | |
|---------------|--------|
| Patient valve | 7 mL |
| Adult mask | 150 mL |
| Child mask | 95 mL |
| Infant mask | 28 mL |

Pressure relief valve

| | |
|----------------------|----------------------------|
| Children and infants | 40 ± 5 cm H ₂ O |
| Adults | 60 ± 5 cm H ₂ O |

| | Bag volume | Stroke volume | Reservoir volume | Suitable body weight |
|--------|------------|---------------|------------------|----------------------|
| Adult | 1600 mL | 700 mL | 2500 mL | > 30 kg |
| Child | 500 mL | 300 mL | 2500 mL | 7-30 kg |
| Infant | 280 mL | 150 mL | 600 mL | < 7 kg |

A stroke volume of 1350 mL can be achieved using two hands.

Minimum cycle rate

| | |
|-------------|-----------------|
| Adult mask | 20 breaths/min. |
| Child mask | 20 breaths/min. |
| Infant mask | 40 breaths/min. |

Oxygen concentration

| | |
|-------------------|--------------------------------------|
| With reservoir | 99% |
| Without reservoir | 45% (models for adults and children) |
| | 90% (model for infants) |

The performance characteristics for manual resuscitators can vary from user to user and are dependent on a variety of factors, such as ambient temperature, patient lung capacity, ventilation frequency and the size of the user's hands.

EN Specifications & Performance Characteristics

Maximum Ventilation Rates

| | Temperature -18°C Ventilation rate | Temperature 22°C Ventilation rate | Temperature 50°C Ventilation rate |
|---------|---------------------------------------|--------------------------------------|--------------------------------------|
| Adult | 20 | 20 | 20 |
| Child A | 30 | 30 | 30 |
| Child B | 20 | 20 | 20 |
| Infant | 60 | 60 | 60 |

The results were obtained under the following conditions:

Adult:

V 600 mL, Expansion capacity 0.02 L/cm H₂O, Resistance 20 cm H₂O/L/s

Child A:

V 70 mL, Expansion capacity 0.01 L/cm H₂O, Resistance 20 cm H₂O/L/s

Child B:

V 300 mL, Expansion capacity 0.01 L/cm H₂O, Resistance 20 cm H₂O/L/s

Infant:

V 20 mL, Expansion capacity 0.001 L/cm H₂O, Resistance 20 cm H₂O/L/s

The correct ventilation frequency may vary. Please follow the current ventilation frequencies as recommended by the AHA.

Delivery Pressure Range

| | |
|--------------------|--|
| Adults: | 60 ± 10 cm H ₂ O maximum (can be deactivated by the user) |
| Children, infants: | 40 ± 5 cm H ₂ O maximum (can be deactivated by the user) |

Stroke Volume Range

| | Using one hand | Using two hands |
|--------|----------------|-----------------|
| Adult | 770 mL | 900 mL |
| Child | 300 mL | 350 mL |
| Infant | 160 mL | 190 mL |

EN Specifications & Performance Characteristics

Oxygen Concentration for Adult Model

The values in brackets refer to ventilation without an oxygen reservoir.
Tidal volume (mL) x ventilation frequency with reservoir (without reservoir).

| Oxygen flow (L/min) | 800 x 12 | 800 x 20 | 850 x 12 | 850 x 20 | 900 x 12 | 900 x 12 |
|---------------------|----------|----------|----------|----------|----------|----------|
| 5 | 47 [31] | 58 [29] | 50 [30] | 63 [29] | 56 [29] | 66 [30] |
| 10 | 71 [35] | 89 [35] | 68 [35] | 87 [35] | 91 [36] | 90 [36] |
| 15 | 71 [35] | 91 [37] | 95 [36] | 95 [36] | 94 [36] | 94 [36] |

Oxygen Concentration for Child Model

The values in brackets refer to ventilation without an oxygen reservoir.
Tidal volume (mL) x ventilation frequency with reservoir (without reservoir).

| Oxygen flow (L/min) | 250 x 20 | 250 x 30 | 300 x 20 | 300 x 30 | 350 x 20 | 350 x 30 |
|---------------------|----------|----------|----------|----------|----------|----------|
| 2 | 64 [34] | 70 [34] | 65 [34] | 62 [33] | 67 [33] | 68 [32] |
| 6 | 98 [45] | 97 [44] | 96 [43] | 96 [41] | 95 [39] | 95 [41] |
| 10 | 98 [46] | 98 [45] | 96 [43] | 97 [44] | 95 [44] | 95 [46] |

Oxygen Concentration for Infant Model


The values in brackets refer to ventilation without an oxygen reservoir.
Tidal volume (mL) x ventilation frequency with reservoir (without reservoir).

| Oxygen flow (L/min) | 160 x 30 | 160 x 40 | 175 x 30 | 175 x 40 | 190 x 30 | 190 x 40 |
|---------------------|----------|----------|----------|----------|----------|----------|
| 2 | 84 [36] | 85 [42] | 82 [36] | 78 [40] | 82 [40] | 72 [40] |
| 6 | 85 [40] | 82 [42] | 82 [36] | 82 [41] | 85 [49] | 74 [44] |
| 10 | 87 [49] | 83 [51] | 85 [36] | 82 [46] | 85 [48] | 77 [52] |

EN The PEEP Valve (Accessory Series)

Terminal Positive Pressure Adjustment Valve (for use in one patient)

- 2 - 10 cm H₂O terminal positive pressure adjustment valve for use in one patient.
- 5 - 20 cm H₂O terminal positive pressure adjustment valve for use in one patient.

Terminal Positive Pressure Adjustment Valve


2 - 10 cm H₂O terminal positive pressure adjustment valve (orange, silicone)

5 - 20 cm H₂O terminal positive pressure adjustment valve (blue, silicone)

Specification

Adjustable range: 2 - 10 cm H₂O and 5 - 20 cm H₂O (+/-2 cm H₂O) Flow adjustment 3 L/min

Materials: PC, silicone, stainless steel

PEEP valve


1. Choose the appropriate PEEP valve within the specification range. (2-10 cm H₂O or 5-20 cm H₂O).
2. Turn the knob of the PEEP valve to the required value indicated on the manometer on the valve base.
3. Attach the PEEP valve to the adapter, squeeze the resuscitator for air exchange and set the required pressure on the PEEP valve.

The PEEP valves are designed for use with the resuscitator to build up positive pressure and expiratory pressure during ventilation. During ventilation, the use of the PEEP valve does not influence the inspiratory resistance or inspiratory oxygen concentration. They can be used during the treatment both to restore spontaneous breathing and for breathing difficulties.

The PEEP valve must be connected to the expiratory flow diverter on the resuscitator.

CE 0123

MD

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Straße 53
50996 Köln · Germany
Tel.: +49 (2236) 3913-0
Fax: +49 (2236) 3913-109
info@dahlhausen.de
www.dahlhausen.de

Rev. 01/2020-03