Technische fiche

Folded FFP2 NR - EN149:2001+A1:2009

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Technische fiche

Dispenser

BruMed
PARTICLE FILTERING
HALF MASK
EN149:2001+A1:2009
Folded FFP2 NR
CE 0598
50
Minimum sales unit: 50 pcs

BruMed
PARTICLE FILTERING
HALF MASK
EN149:2001+A1:2009
Folded FFP2 NR
CE 0598
50 PCS
Minimum sales unit: 50 pcs

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EN Particulate Respirator. See instructions for use. For use against solid and liquid aerosols only. For single use only.

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| Marking on Product | Description on label | Explanation |
|---------------------|----------------------|-----------------------------|
| | | Identification Mask |
| Folded FFP2 NR | Folded | Product identification |
| EN 149:2001+A1:2009 | CE 0598 | CE mark |
| | EN149:2001+A1:2009 | Number of European standard |
| | FFP2 NR | Protection category |

Eu -Type Examination Notified Body
Notified Body: SGS Fimko Oy
Address: Takomotie 8 FI-00380 Helsinki, Finland.
Notified Body No: 0598

Information of Manufacturer
Manufacturer: Sample King Manufacturing Ltd
The EU declaration of conformity accompany with product

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Usage / Limitations
Filtering half masks for single use. NR = not reusable for a maximum of 1 working day of 8 hours.
This respirator is suitable for use in protection against non-toxic solid and liquid aerosols. Do not use out of the scope or use defined in the warnings. Failure to properly use this product may result in serious health damage or death.
- FFP1 NR: Filter Efficiency 80%; Allocated Protection Factor (FPA) is 4; Examples of applications are handling of stone / rubble / cellulose.
- FFP2 NR: Filter Efficiency 94%; Allocated Protection Factor (FPA) is 10; Examples of applications are sanding of soft wood, composite materials, rust plaster, plastics / cutting, deburring, grinding, drilling of metal.
- FFP3 NR: Filter Efficiency 99%; Allocated Protection Factor (FPA) is 20; Examples of applications are sanding of hard wood (beech, oak) / treatment of wood using copper, chrome or arsenic based products / impact stripping of paint / sanding of cement.