
SINCE 1966  
MADE IN ITALY

# Air Pro

by **FLAEM**®

**3000**  
Plus Version

- IT > ISTRUZIONI PER L'USO** - Aerosol professionale per uso ospedaliero, home care e noleggio  
**EN > INSTRUCTIONS FOR USE MANUAL** - Professional aerosol therapy unit for hospital use, home care and renting  
**FR > MODE D'EMPLOI** - Générateur d'aérosol professionnel à usage hospitalier, de soin à domicile et de location  
**NL > GEBRUIKSAANWIJZING** - Professionele aërosol voor gebruik in ziekenhuizen, home care en verhuur  
**DE > BEDIENUNGSANLEITUNG** - Professioneller Inhalator für den Gebrauch im Krankenhaus, Homecare und zum Mieten  
**ES > MANUAL DE INSTRUCCIONES** - Nebulizador profesional para uso hospitalario, asistencia domiciliaria y alquiler  
**PL > INSTRUKCJA UŻYCIA** - Profesjonalny aerzol do użytku w szpitalach, domach opieki i wypożyczalniach  
**EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ** - Επαγγελματικό αερόλυμα για νοσοκομειακή χρήση, φροντίδα στο σπίτι και ενοικίαση

**SCHEMA DI COLLEGAMENTO PER NOMENCLATURA - CONNECTION DIAGRAM FOR NOMENCLATURE**  
**SCHEMA DE CONNEXION POUR NOMENCLATURE - NOMENCLATUUR VERBINDINGSSCHEMA**  
**SCHALTPLAN FÜR DIE BEZEICHNUNGEN - ESQUEMA DE CONEXIÓN PARA NOMENCLATURA**  
**СХЕМА СОЕДИНЕНИЯ СОСТАВНЫХ ЧАСТЕЙ - SCHEMAT POŁĄCZEŃ DO NOMENKLATURY -**  
**ΣΧΕΔΙΑ ΣΥΝΔΕΣΗΣ ΓΙΑ ΟΝΟΜΑΤΟΛΟΓΙΑ**

ITALIANO  
Pag. 1

ENGLISH  
Pag. 10

FRANÇAIS  
Pag. 19

NEDERLANDS  
Pag. 28

DEUTSCH  
Pag. 37

ESPAÑOL  
Pag. 46

POLSKI  
Pag. 55

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## Mod. P0712EM F3000 Plus Version

### Aerosol therapy apparatus

We are pleased you have purchased our product and we thank you for your trust in us.

Read these instructions carefully and retain them for future reference.

Only use the apparatus as described in this manual.

## INTENDED USE

**Air Pro 3000 Plus Version** is a professional medical device whose intended use is the administration of medication prescribed or recommended by your doctor upon assessing the patient's general conditions, for the treatment of breathing ailments. It can be used in health care structures like hospitals, in home care and for hire. The product should be used by legally authorized medical personnel (doctors, nurses and therapists), by the home therapist, or by the patients. It is highly recommended that the health professional and/or the patient read and understand the use and maintenance information of this manual.

Carefully follow the user and maintenance instructions to ensure proper functioning and to extend the life of the device.

**Air Pro 3000 Plus Version** is fitted with a knob to adjust the pressure (A3), and with pressure gauge (A6) to display the set pressure value.

The device does not require lubrication, is manageable, simple to use, reliable, resistant and quiet.

Note: Use only original FLAEM accessories.

## STANDARD EQUIPMENT INCLUDES:

### A Aerosol Device (main unit)

- A1- Switch
- A2- Air outlet port
- A3- Pressure adjustment knob
- A4- Nebulizer port
- A5- Carrying handle
- A6- Pressure gauge
- A7- Power supply cord socket
- A8- Power cord

### B Connection tube (main unit / nebuliser)

### C Accessories

- C1- RF7 Dual Speed Plus Nebulizer
  - C1.1- Lower part
  - C1.2- Nozzle
  - C1.3- Upper part
  - C1.4- Speed selector with valve system
- C2 - Mouthpiece with valve
  - C2.1- Expiratory valve
- C3- Non-invasive adult nasal prong
- C4- Child mask
- C5- Adult mask
- C6- Nebulization manual control

## ⚠ IMPORTANT NOTICES

- This device is also intended for direct use by the patient.
- Before using the product for the first time, and periodically during its lifetime, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised FLAEM service centre or to your trusted dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- The expected medical life of the accessories is 1 year. It is, however, advisable to replace the nebuliser cup every 6 months in the event of intense use (or earlier if the cup is obstructed) to always guarantee maximum therapeutic efficacy.
- Children and people who are not self-sufficient must use the device under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.
- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Always keep the power supply cable away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the device with wet hands. Do not use the device in damp environments (for example, while taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.
- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.

- Do not expose the device to particularly extreme temperatures.
- Do not place the device near sources of heat, in direct sunlight or in excessively hot rooms.
- Do not obstruct or put objects into the filter or its related housing in the device.
- Never obstruct the air vents located on both sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Make sure there is no material obstructing the air vents before each use.
- Do not put any objects in the air vents.
- Repairs, including the replacement of the supply cord, are to be carried out by FLAEM authorised personnel only, by complying with the information provided by the manufacturer.
- The average operating life of the compressor is 2000 hours approximately.
- WARNING: Do not modify this device without authorisation from the manufacturer.
- The manufacturer, the vendor and the importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- Interactions: the materials used in contact with medication have been tested with a vast range of medications. However, in view of the variety and continuous evolution of pharmaceuticals, interactions cannot be ruled out. We recommend using the medication as soon as possible once it has been opened and preventing prolonged exposure in the nebuliser cup. The materials used in the medical device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- The amount of time required between preservation conditions and use is approximately 2 hours.

## INSTRUCTIONS FOR USE


**Before each use, clean hands thoroughly and clean the device as described in the section on “CLEANING, SANITISATION, DISINFECTION, STERILISATION”. During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination.**

**This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.**


1. Insert power supply cord (A8) into socket (A7) and then connect it to a mains outlet adequate to the device voltage. The position of the socket must be such that the device can be easily unplugged from the mains network



2. Pour the medication prescribed by the doctor into the bottom part (C1.1)

3. Close the nebuliser by turning the top (C1.3) clockwise.

4. Connect one end of the connection pipe (B) to the device. Connect the other end of the pipe (B) to the manual nebulisation control (C6). Connect manual nebulisation control (C6) to the atomizer (C1).

5. Connect accessories as indicated in the "Connection diagram" on the cover. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

## PRESSURE ADJUSTMENT

### 6.a Operating pressure

After connecting nebuliser RF7 Plus to the device with the Selector Button (C1.4) closed in the position With Valve, operating pressure can be adjusted by rotating Knob (A3), clockwise to increase it or anti-clockwise to decrease it, until selecting the colour sector visible on the Pressure Gauge (A6) matching the desired particle size shown in the table below.

Throughout the device life, from time to time it is necessary to check maximum pressure attainable with the RF7 Plus nebuliser connected and with Knob (A3) turned on MAX.

If the maximum value of the green sector in the Pressure Gauge is not reached (A6), the device must be serviced.

## CHOICE OF DELIVERY METHOD

### 6.b Delivery method


Nebulisation can be obtained by maintaining the valve system active, or selecting the delivery mode at top speed, by means of the Selector Button (C1.4) of RF7 Plus atomiser. The mode must be selected after setting the operating pressure.

7. Start the device by means of switch (A1) and breathe deeply in and out.

8. Upon completing application, turn off the device and unplug it.


9. If after the therapy session an evident deposit of moisture forms within the pipe (B), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.

## RF7 PLUS TECHNICAL NEBULISATION FEATURES


	Green Area	Yellow Area	Orange Area	Green Area	Yellow Area	Orange Area
Pressure gauge value bar	1,15 ÷ 1,45	0,90 ÷ 1,15	0,70 ÷ 0,90	1,15 ÷ 1,45	0,90 ÷ 1,15	0,70 ÷ 0,90
Nozzle flow (1) l/min'	8.2	7.1	6.8	8.2	7.1	6.8
Selector Button Position	with Valve	with Valve	with Valve	MAX	MAX	MAX
MMAD $\mu\text{m}$ (2)	2.21	2.76	2.95	2.38	2.91	2.44
(2) Breathable fraction < 5 $\mu\text{m}$ (2)	84.2%	78.5%	74.7%	81.5%	76.0%	80.7%
Delivery ml/min' (1) approx	0.29	0.23	0.18	0.65	0.53	0.42
(1) Data recorded according to internal Flaem procedure I29-P07.5. Delivery Speed values are subject to variations depending on the patient's respiratory ability. (2) In vitro testing performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma and in compliance with the EN 13544-1: 2007 + A1 European Standard for aerosol therapy devices Standard. More details are available on request.						


## SOFT TOUCH MASKS

The new SoftTouch Masks feature an adaptive edge made of a soft biocompatible material ensuring an excellent adherence to the face and a Dispersion Limiting Vent. These innovations, which make them stand out from the rest, allow increased sedimentation of medication in the lungs and reduced dispersion into the air.


During the inspiratory phase, the tab, acting as a Dispersion Limiting Device, bends inwards towards the mask


During expiration the vent bends outwards.

## USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (C6), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (C6) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

## CLEANING, SANITISATION, DISINFECTION, STERILISATION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

### DEVICE (A) AND TUBING EXTERIOR (B)

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

### ACCESSORIES

Open the nebuliser by turning the upper part (C1.3) anticlockwise, remove the nozzle (C1.2) and the speed selector with valve system (C1.4) from the upper part (C1.3) as shown in the "Assembly diagram" in point C1. Then proceed according to the following instructions

### CLEANING AT HOME - SANITISATION AND DISINFECTION

#### SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

**method A:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

**method B:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** in the dishwasher with a hot cycle.

**method C:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40°C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

#### DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

**method A:** Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

#### Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.

- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

**method B:** Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

**method C:** Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

### CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

**method A:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

**method B:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** in the dishwasher with a hot cycle.

#### DISINFECTION

Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

#### Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

If you want to also perform the STERILISATION, jump to the STERILISATION paragraph.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

#### STERILISATION

Accessories **C1 (Assembled nebulizer)-C2-C3-C4-C5** can be sterilised.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

**Implementation:** Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and a time of 10 minutes first.

**Storage:** Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging. The sterilisation procedure is validated in its conformity to ISO 17665-1.


**At the end of each use store the device complete with accessories in a dry place away from dust.**




### AIR FILTERING

The device is fitted with a suction filter, located inside the pressure-adjustment knob, which must be replaced when it is dirty or discoloured. To replace it, lift the pressure-adjustment Knob (A3) and extract the Filter.

Do not wash or reuse the same Filter. Regular filter replacement is necessary to assure proper compressor performance. The Filter must be regularly checked. The filter was designed to always remain secured in its seat. Do not replace the filter during use.

**Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories.**

SPARE PARTS AND ACCESSORIES	DESCRIPTION	CODE	PACKING
	RF7 Dual Speed plus nebulizer + mouthpiece with valve + nose piece	ACO464P	(1PC)
	RF7 dual speed plus nebulizer	ACO469P	(50PCS)
	2m pipe with fittings	ACO578P	(1PC)
	2m pipe with fittings	ACO735P	(50PCS)
	Mouthpiece with valve	ACO441P	(50PCS)
	Non-invasive nose piece	ACO448P	(50PCS)
	Adult mask	ACO473P	(50PCS)
	Paediatric mask	ACO472P	(50PCS)
	Complete newborn accessory set	ACO392P	(1PC)
	Newborn accessory set size 2	ACO391P	(1PC)
	Newborn accessory set size 1	ACO390P	(1PC)
	Newborn accessory set size 0	ACO389P	(1PC)
	PP adult nose piece	ACO450P	(50PCS)
	PP paediatric nose piece	ACO451P	(50PCS)
	Porex filter	ACO45	(50PCS)
	RF6P nebulizer kit w/anti-dispersal filter	ACO388P	(1PC)
	Single-use filter for anti-leakage kit	ACO396P	(50PCS)

## TROUBLE-SHOOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power cable has not been correctly inserted into the socket of the device or the power socket	Correctly insert the power cable in sockets.
The device does not nebulise or nebulises insufficiently.	The medication has not been inserted in the nebuliser.	Pour the right amount of medication in the nebuliser.
	The nebuliser has not been properly fitted.	Disassemble the nebuliser and reassemble it correctly as shown in the connection diagram on the cover.
	The nebuliser nozzle is clogged.	Disassemble the nebuliser, remove the nozzle and perform cleaning operations Failure to remove medication deposits from the nebuliser affects its efficiency and operation Strictly comply with the instructions contained in the CLEANING, SANITISATION AND DISINFECTATION chapter
	The air pipe is not correctly connected to the device	Ensure the device's accessories are properly connected to the air intake (see connection diagram on the cover).
	The air pipe is bent, damaged or twisted	Unwind the pipe and ensure it is not crushed or punctured. If necessary replace it.
	The air filter is dirty	Replace the filter
The apparatus is noisier than usual	The pressure regulator is not correctly inserted in its housing	Fully and properly insert the pressure regulator in its housing

If after verifying the above mentioned conditions the device should not operate properly, contact your trusted retailer or the nearest authorised FLAEM service centre.

## SYMBOLS

### Certification TÜV

 CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments

 Class II device

 Important: check the operating instructions

 "OFF" for part of equipment

 "ON" for part of equipment

 Complies with: European standard EN 10993-1 "Biological Evaluation of Medical Devices" and European Directive 93/42/EEC "Medical Devices". Phthalates free. In conformity with Reg. (EC) no. 1907/2006

 Temperature limit

 Atmospheric pressure limit

 Serial number of device

 Manufacturer

 Type BF applied part

 Alternating current

**IP21** Enclosure protection rating: IP21.  
(Protected against solid bodies over 12 mm. Protected against access with a finger.  
Protected against vertically falling water drops.)

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### APPLIED PARTS

The applied parts of type BF are: (C2, C3, C4, C5)

### TECHNICAL SPECIFICATIONS RF7 Dual Speed Plus Nebulizer

Minimum fill volume:	2 ml
Maximum fill volume:	8 ml

## ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website [www.flaemnuova.it](http://www.flaemnuova.it). The appliance may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnoses or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

## DEVICE DISPOSAL

 In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

## PACKAGING DISPOSAL

Product box and accessories box

Product packaging bag and accessories packaging bag

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© 2021 FLAEM NUOVA®
**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**
FLAEM NUOVA S.p.A.