# CORE STABILIZERS

## Pressure Bio-Feedback Unit

## Operating Instructions

*ENGLISH.....page 3*

*DEUTSCH.....page 15*

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*ESPAÑOL.....page 39*

*ITALIANO.....page 52*

*NEDERLANDS.....page 63*

*DANSK.....page 75*

*SVENSKA.....page 87*

# TABLE OF CONTENTS

ENGLISH .....	3
DEUTSCH .....	15
FRANÇAIS .....	27
ESPAÑOL .....	39
ITALIANO .....	52
NEDERLANDS .....	63
DANSK .....	75
SVENSKA .....	87

The Pressure Bio-Feedback Unit (PBU), designed by physical therapists, is a simple device which registers changing pressure in an air filled pressure cell. This allows body movement, especially spinal movement, to be detected during exercise. The unit consists of a combined gauge/inflation bulb connected to a pressurecell.

Core Stabilization training involves the re-learning of coactivation of co-contraction patterns of the transverse abdominis and the lumbar multifidus to provide local spinalsegmental support. Published research has shown that when functioning normally, the transverse abdominis and the lumbar multifidus act in concert, increasing tension on the thoracolumbar fascia acting like a corset, providing stability to the lumbar spine.

Research has shown these types of exercise are especially important for the prevention and treatment of low back and neck pain (of various pathologies). The exercise techniques discussed in this instruction manual involve slow and controlled body movement. The PBU is used to monitor and provide feedback on body movement during exercise. This promotes effective exercise to improve back and neck pain.

The concepts put forth in this literature were in effect at the time of publication. However, due to the company's policy of continuous improvement, changes may be made at any time without obligation on the part of the designing physical therapists or the company.

## WARNING

**CAUTION:** This product contains **natural rubber latex** which may cause allergic reactions.

# CORE STABILIZERS

## Pressure Bio-Feedback Unit

**A simple device that provides feedback to ensure quality and precision in exercise performance and testing.**

The PBU is useful in an infinite number of exercises for all parts of the body. This instruction manual contains information on its use for exercises aimed at enhancing the stability of the lumbo-pelvic and cervical region.

Considerable research has been completed on the specific stabilization exercises related to these areas.

### GENERAL INSTRUCTIONS

- The three-chamber pressure cell of the Core stabilizers is placed between the part of the body requiring monitoring and a firm surface (e.g., floor, back of chair, plinth, bed, wall).
- Tighten the screw (gently) at the base of the pressure gauge. Pump to inflate the three chambers of pressure cell until it molds between the body part and the supporting surface. A pressure of 40 mmHg (orange band) is suitable for the resting pressure of the inflated cell.
- The initial baseline pressure may decrease initially as the air backflows into the tube. Reinflate before commencing exercise.
- Changes in body weight on the cell on any of the three compartments will register a pressure change on the gauge.
- Pressure on the cell can be increased, decreased or maintained depending on the exercise/movement required.
- After exercise, air is released by loosening the screw.

# PRECISION IN SPECIFIC MUSCLE TESTING

The deep muscle system in the lumbar spine and pelvic region is directly responsible for stabilizing the vertebral segments and sacro-iliac joints. The larger, more superficial muscles are involved in moving the trunk. The deep muscles are underactive in low back pain and the more superficial muscles are often overactive. In order to retrain the deep muscles, specific exercise and testing with the PBU is required.

## 1. The Prone Test for Transversus Abdominis and Internal Oblique

- Place the three-chamber pressure cell under the abdomen and inflate to baseline of 70 mmHg. (brown band)
- Draw abdominal wall up and in without moving the spine or pelvis.
- Pressure should decrease 6-10 mmHg.
- Hold 10-15 seconds, breathe normally.
- Perform 10 repetitions.

Lying on abdomen

## 2. Training the Corset Action of Transversus Abdominis in Supine

- Place the three-chamber pressure cell under the lumbar spine and inflate to a baseline of 40 mmHg (orange band).
- Draw in the abdominal wall without moving the spine or pelvis.
- Pressure should remain at 40 mmHg (i.e. no movement of the spine).
- Hold for 10-15 seconds; breathe normally.
- Perform 10 repetitions.

Lying on back

## 3. Training the Corset Action of Transversus Abdominis with Leg Loading

- Place the three-chamber pressure cell behind the lumbar spine and inflate to baseline of 40 mmHg (orange band).
- Draw in the abdominal wall without moving the spine or pelvis.
- Pressure should remain at 40 mmHg (i.e. no movement of the spine) while leg lifts.
- Hold for 10-15 seconds; breathe normally.
- Repeat 10 times with each leg.

Standing

Lying (controlled leg movement)

## 4. Lumbar Stabilization During Stretching Techniques

The PBU helps to stabilize the spine while adjacent body segments are moved in order to stretch the appropriate muscles.

### To Test/Stretch:

- Hip Flexors (Thomas Test)
- Rectus Femoris
- Latissimus Dorsi and Pectoralis Major.

Follow stabilization procedure and maintain pressure during the stretch.

Test/stretch for tight hip flexors  
(Thomas Test)

Test/stretch for tight  
Rectus Femoris

Test/stretch for tight  
Latissimus Dorsi and  
Pectoralis Major

### To Test/Stretch: Tensor Fascia lata

- Have patient lying on their side.
- Place pressure cell between lateral trunk (lumbar spine level) and supporting surface.
- Inflate pressure cell to 40 mmHg. (orange band)
- Maintain lumbo-pelvic stability by maintaining pressure at 40 mmHg during the stretching maneuvers.

Test/stretch for tight  
tensor fascia lata

## 5. Lower Trapezius

- Place the three-chamber pressure cell under abdomen and inflate to baseline of 70 mmHg. (brown band)
- Draw abdominal wall up and in as described in exercise 1. (Prone test for Transverse Abdominis and Internal Obliques.)
- Lift one arm and draw the scapula down and in, towards the spine. The pressure must remain **constant**.
- Hold for 5 seconds.
- Perform 10 repetitions.


**Lying on abdomen  
(arm exercise)**

**The PBU is used to monitor proximal trunk stability to detect substitution for the tested muscle's action.**

## **6. Iliopsoas**

- Place the three-chamber pressure cell behind the lumbar spine and inflate to baseline of 40 mmHg. (orange band)
- Draw in abdominal wall without thoracic or pelvic movement.
- Pressure should increase 8-10 mmHg.
- Lift one leg slowly to test position with the other leg unsupported (i.e., no stool). Pressure should remain constant at 48-50 mmHg.
- Hold for 5 seconds.
- Perform 10 repetitions.

## **7. Gluteus Maximus**

- Place the three-chamber pressure cell under the abdomen and inflate to baseline of 70 mmHg. (brown band)
- Draw abdominal wall up and in as described in Exercise 1. (Prone test for Transverse Abdominis and Internal obliques.)
- Extend hip by lifting knee up 5 cm (2 in) off supporting surface. The pressure must remain constant.
- Hold for 5 seconds.
- Perform 10 repetitions.

## 8. Training of the Deep Cervical Flexors

- Position the folded cell (secure with pressure studs) under the neck so that it abuts against the occiput. Do not let it slide down to the lower cervical area.
- Inflate to a baseline of 20 mmHg. (red band)

NOTE: Do not inflate the PBU before insertion behind the neck. Because the air has to completely enter all of the bag's pressure cells, it is helpful to squeeze the pressure bag to speed up the air distribution. Reinflate and repeat the gentle squeeze. This may take two to three repetitions to complete the process.

- The movement the patient is to perform is a gentle nodding of the head as if they were saying "yes".
- Instruct the patient to place the tongue on the roof of the mouth, lips together but teeth just separated. This will discourage substitution with the platysma or hyoids.
- Instruct the patient to gently nod to target 22 mmHg, just one mark on the pressure dial. (green band) See if the patient can hold the position steadily. If successful, relax and repeat at each target position of 24 mmHg through 30 mmHg. (yellow, blue and gray bands)
- The pressure that the patient can hold steady, with minimal superficial muscle activity, is the one on which you will measure endurance capacity. (i.e. 10 repetitions held for 10 seconds each)

Lying on back (neck exercise)

# MONITORING OF LUMBAR SPINE STABILIZATION FOR SAFETY AND PRECISION OF EXERCISE AND PROTECTION OF LUMBAR SPINE

## **Limb Loading**

- PBU is used to monitor the appropriate use of the muscles required in the stabilization and protection of the lumbar spine during lower and upper limb testing and exercise.

## **Stabilization Procedure**

- Place the three-chamber pressure cell longitudinally under the back of the side on the spine opposite the limb being lifted in exercise and inflate to baseline of 40 mmHg. (orange band)
- Control the back position during limb loading activities and keep the pressure constant.
- In addition, a presetting action of the abdominals may be performed prior to limb loading.

## **Either**

Presetting as for Exercise 2. The pressure will increase 8-10 mmHg. Keep pressure constant.

## **OR**

Presetting with a pelvic tilt position. The pressure will increase 20-50 mmHg. Keep pressure constant.

Lying (controlled leg movement)

Standing

## SPECIFICATIONS

<b>MEASURING RANGE</b>	Analog 0-200 mmHg pressure
<b>ACCURACY</b>	$\pm 3$ mmHg pressure
<b>INFLATE</b>	Squeeze bulbs (after tightening screw valve)
<b>DEFLATE</b>	Loosen screw valve

## RECOMMENDATIONS

1. The connector tube should be disconnected from the dial gauge/inflation bulb before storage.
2. The three-chamber pressure cell must be stored flat or folded along the seams when not in use.
3. Cleaning the PBU should be done only with a damp cloth.

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