DE Diagnosewaage Gebrauchsanweisung	2	IT Bilancia diagnostica Istruzioni per l'uso	37
EN Diagnostic bathroom scale Instruction for Use.....	11	TR Diyagnoz terazisi Kullanma Talimatı	46
FR Pèse-personne impédancemètre Mode d'emploi	19	RU Диагностические весы Инструкция по применению ...	54
ES Báscula diagnóstica Instrucciones para el uso	28	PL Waga diagnostyczna Instrukcja obsługi	63

Inhalt

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Germany

Contents

1. Getting to know your instrument 11	8. Taking measurements..... 14
2. Signs and symbols 11	9. Evaluation of results 15
3. Safety notes..... 12	10. Replacing batteries..... 17
4. Unit description 13	11. Cleaning and care of the unit 17
5. Initial use 13	12. Disposal..... 17
6. Information 13	13. What if there are problems? 18
7. Settings 14	14. Warranty/service..... 18

Dear Customer

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain. This diagnostic scale is designed to make a significant contribution to your health.

With kind regards
Your Beurer team

1. Getting to know your instrument

Function of the unit

The purpose of this digital scale is weighing and calculating your personal fitness data.

It is intended for private use.

The scale is equipped with the following functions, which can be used by up to 10 different people:

- body weight measurement,
- measurement of body fat percentage,
- body water percentage,
- muscle percentage,
- bone mass,
- basic and active metabolic rate.

This scale also has the following functions:

- interpretation of body fat values,
- switch between kilograms “kg”, pounds “lb” and stones “st”,
- automatic shutoff function,
- battery change indicator for weak batteries.

2. Signs and symbols

The following symbols appear in these instructions.

Warning

Warning instruction indicating a risk of injury or damage to health.

Important

Safety note indicating possible damage to the unit/accessory.

Note

Note on important information.

3. Safety notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

WARNING

- **The scales must not be used by persons with medical implants (e.g. heart pacemakers) Otherwise their function could be impaired.**

- Do not use during pregnancy.
- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).
- Caution! Do not step onto the scale with wet feet or when the scale's surface is damp – danger of slipping!

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
- **⚠ Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
- **⚠ Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.


General notes

- The unit is for personal use only and is not intended for medical or commercial applications.
- Note that technical tolerances mean that results may vary, because the scale is not calibrated for professional medical use.
- The capacity of the scale is 150 kg (330 lb, 23 st:8 lb). For weight measurement and bone mass measurement, the results are shown in 100 g increments (0.2 lb).
- The measuring results of the body fat, body water and muscle percentage are shown in 0.1 % increments.
- The caloric requirements is specified in steps of 1 kcal.
- On delivery, the scale is set to "cm" and "kg". At the back of the scales is a switch to change between "pounds" and "stones".
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- The scale should be cleaned occasionally with a damp cloth. Do not use abrasive detergents and never immerse the scale in water.
- Protect the unit against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stoves, heating radiators).
- Repairs may only be carried out by Beurer customer service or authorised dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- Please contact the specified service address to obtain further information, such as the CE Declaration of Conformity.
- Please note that the sensor buttons react to conductive objects such as metal or water-filled containers etc.,

4. Unit description

Overview

1. Electrodes
2. Display
3. “Up” key
4. “SET” key
5. “Down” key

5. Initial use

Inserting batteries

If provided, pull out the battery insulating strip at the battery compartment or remove the protective film from the battery and insert the battery with the correct polarity. If the scale fails to operate, remove the battery completely and reinsert it.

Changing the weight unit

On delivery, the scale is set to “cm” and “kg”. At the back of the scales is a switch with which you can switch between “inches”, “pounds” and “stones” (lb, st).

Setting up the scale

Place the scale on a firm level floor (no carpet); a firm floor covering is required for correct measurement.

Switch on the scale

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor.

6. Information

The measuring principle

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, gender and degree of activity.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

General tips

- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist.
Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.
- Stand still and upright during the measurement.
- Wait a few hours after unaccustomed physical activity.
- Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

Limitations

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women,
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

7. Settings

Enter your personal information before using the scale.

On delivery, the scale is set to “cm” and “kg”. At the back of the scales is a switch with which you can switch between “inches”, “pounds” and “stones” (lb, st).

Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 10 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later.

- Switch on the scale. Wait until the display shows “0.0 kg”.
- Start setting with the SET key. The first memory space flashes in the display.
- Select the required memory space with the ▲/▼ key and confirm with the SET key.
- Now you can adjust the following settings:

User data	Parameters
Memory position	1 to 10
Body height	100 to 220 cm (3' 3.5" to 7' 2.5")
Age	10 to 100 Jahre
Sex	man (♂), woman (♀)
Degree of activity	1 to 5

- Modifying values: Press ▲ or ▼ or press and hold for more rapid progress.
- Confirming entries: Press the SET button.
- Once the values have been displayed, “0.0 kg” appears in the display.
- Then the scale is ready for measuring. If you do not carry out a measurement, the scale switches off automatically after several seconds.

Degrees of activity

Selection of the degree of activity must refer to the medium and long term.

Degree of activity	Physical activity
1	None.
2	Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises).
3	Medium: Physical effort for 30 minutes at least 2 to 4 times a week.
4	High: Physical effort for 30 minutes at least 4 to 6 times a week.
5	Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily.

All parameters must be entered before the body fat and other values can be entered.

8. Taking measurements

Place the scale on a firm level floor (no carpet); a firm floor covering is required for correct measurement.

Only weight measurement

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor.

“8888” (Fig. 1) will appear on the display as a self test until “0.0 kg” is displayed (Fig. 2). Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight.

Fig. 1

Fig. 2

Weighing, carrying out diagnostics

- Switch on your scale and wait until “0.0 kg” appears (Fig. 2).
- Select the memory space where your personal data are stored by repeatedly pressing the ▲ or ▼ keys. These are then displayed in sequence until “0.0 kg” appears.
- Step onto the scale with bare feet and make sure that you are standing still with equal weight distribution. Stand on the stainless steel electrodes with both feet.

Note: Your feet, legs, calves and thighs must not touch each other.

Otherwise the measurement cannot be performed correctly. The scale immediately begins to measure your weight.

The scales begin weighing immediately. First, the weight is displayed.

Whilst the further parameters are being measured an “0000” appears (Fig. 3).

Your measured weight is displayed shortly after this.

Fig. 3

The following appears:

1. Weight in kg
2. Body fat in % with interpretation of BF
3. Body water in % ≈
4. Muscle percentage in % ⇄⇄
5. Bone mass in kg ⇄
6. Basal metabolic rate in kcal
7. Active metabolic rate in kcal

Switching off the scale

The scale automatically switches off after approx. 10 seconds.

9. Evaluation of results

Body fat percentage

Body fat is evaluated and displayed by a bar below the body fat display.

The following body fat percentages are guide values (contact your physician for further information).

Man

Age	low	normal	high	very high
10–14	<11 %	11–16 %	16.1–21 %	>21.1 %
15–19	<12 %	12–17 %	17.1–22 %	>22.1 %
20–29	<13 %	13–18 %	18.1–23 %	>23.1 %
30–39	<14 %	14–19 %	19.1–24 %	>24.1 %
40–49	<15 %	15–20 %	20.1–25 %	>25.1 %
50–59	<16 %	16–21 %	21.1–26 %	>26.1 %
60–69	<17 %	17–22 %	22.1–27 %	>27.1 %
70–100	<18 %	18–23 %	23.1–28 %	>28.1 %

Woman

Age	low	normal	high	very high
10–14	<16 %	16–21 %	21.1–26 %	>26.1 %
15–19	<17 %	17–22 %	22.1–27 %	>27.1 %
20–29	<18 %	18–23 %	23.1–28 %	>28.1 %
30–39	<19 %	19–24 %	24.1–29 %	>29.1 %
40–49	<20 %	20–25 %	25.1–30 %	>30.1 %
50–59	<21 %	21–26 %	26.1–31 %	>31.1 %
60–69	<22 %	22–27 %	27.1–32 %	>32.1 %
70–100	<23 %	23–28 %	28.1–33 %	>33.1 %

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10–100	<50 %	50–65 %	>65 %

Woman

Age	poor	good	very good
10–100	<45 %	45–60 %	>60 %

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

Age	low	normal	high
10–14	<44 %	44–57 %	>57 %
15–19	<43 %	43–56 %	>56 %
20–29	<42 %	42–54 %	>54 %
30–39	<41 %	41–52 %	>52 %
40–49	<40 %	40–50 %	>50 %
50–59	<39 %	39–48 %	>48 %
60–69	<38 %	38–47 %	>47 %
70–100	<37 %	37–46 %	>46 %

Woman

Age	low	normal	high
10–14	<36 %	36–43 %	>43 %
15–19	<35 %	35–41 %	>41 %
20–29	<34 %	34–39 %	>39 %
30–39	<33 %	33–38 %	>38 %
40–49	<31 %	31–36 %	>36 %
50–59	<29 %	29–34 %	>34 %
60–69	<28 %	28–33 %	>33 %
70–100	<27 %	27–32 %	>32 %

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise.

With appropriate muscle building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water). Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender). no recognised guidelines or recommendations relating to bone mass measurement.

ATTENTION:

Please do not confuse bone mass with bone density.

Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest to maintain its basic functions (e.g. while lying in bed for 24 hours). This value largely depends on weight, height and age. It is displayed on the diagnostic scale in kcal/day units using the scientifically recognized Harris-Benedict formula.

Your body requires this amount of energy in any case and it must be reintroduced into your body in the form of nutrition. If you take on less energy over the longer term, this can be harmful to your health.

AMR

The active metabolic rate (AMR) is the amount of energy required daily by the body in its active state. The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1–5).

To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will decrease. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

-  Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

The interpretation of the results will depend on changes in your overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Sudden changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass.
- Ideally you should support your diet with physical activity, fitness or power training. This way you can increase your muscle percentage over the medium term.
- Body fat, body water or muscle percentages should not be added (certain elements of muscle tissue also contain body water).

10. Replacing batteries

Your scale is provided with a “Change battery” indicator. If the scale is operating with a battery which is too weak, “Lo” appears on the display and the scale will switch off automatically. In this case, the battery must be exchanged (3 x 1.5 V AAA).

 **NOTE:**

- When changing batteries, always use batteries of the same type, brand and capacity.
- Do not use rechargeable batteries.
- Use batteries that are free of heavy metals.

11. Cleaning and care of the unit

The unit should be cleaned occasionally.

Clean using a damp cloth, to which you can apply a little detergent if necessary.

IMPORTANT

- Do not use aggressive solvents or cleaning agents!
- Never immerse the unit in water!
- Do not wash the unit in a dishwasher!

12. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner 

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,
 Cd = Battery contains cadmium,
 Hg = Battery contains mercury.

13. What if there are problems?

If the scale detects an error when measuring, the following is displayed.

Display	Cause	Remedy
Err 1	The maximum load-bearing capacity of 150 kg was exceeded.	The load must not exceed 150 kg (330 lb, 23 st : 8 lb)
	The fat percentage is outside the measurable range (less than 5 % or greater than 50 %).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.
Lo	The scale batteries are almost flat.	Replace the batteries of scale.
Err 2	Scale zero setting is incorrect.	Wait until the scale switches itself off. Now correctly activate scale, wait for "0.0 kg" to appear and repeat measurement.
8888	Self test.	Place scale on an even, hard surface and wait for "0.0 kg" to appear.

14. Warranty/service

In case of a claim under the warranty please contact your local dealer or the local representation which is mentioned in the list "service international".

In case of returning the unit please add a copy of your receipt and a short report of the defect.

The following warranty terms shall apply:

1. The warranty period for BEURER products is either 5 years or- if longer- the country specific warranty period from date of Purchase.
 In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Repairs (complete unit or parts of the unit) do not extend the warranty period.
3. The warranty shall not be valid for damages because of
 - a. improper treatment, e.g. nonobservance of the user instructions.
 - b. repairs or tampering by the customer or unauthorised third parties.
 - c. transport from the manufacturer to the consumer or during transport to the service centre.
 - d. The warranty shall not be valid for accessories which are subject to normal wear and tear (cuff, batteries etc.).
4. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

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