


**HOSPITAL CARE**

# Surshield™ Versatus -W

Winged I.V. Catheter with Protector  
**Product Information**



The Surshield protector allows to completely close the needle, in order to prevent accidental needle stick injuries when disposing of the intravenous cannula. Using a safety catheter does not require any change in the insertion technique.

Description		Raw Materials
1	Catheter	Polyurethane + Ba SO <sub>4</sub>
2	Catheter hub Caulking pin	Polypropylene + master batch Polyoxymethylene
3	Cannula (internal needle)	Stainless steel
4	Needle hub	Acrylonitrile butadiene styrene
5	Luer-lock cap	Polypropylene
6	Vent fitting filter Vent fitting body	Polyethylene Acrylonitrile butadiene styrene
	Protector	Polyethylene
Surshield Protector		
7	Protection housing	Polyamide
	Safety fingers	Polysulfone
	Sliding bush + Needle blocker	Polycarbonate + stainless steel cup

Does not contain components made of natural rubber latex/PVC.

Designed and  
 manufactured by skilled  
 technicians in Europe

## Your catheters of choice, part of a wider I.V. portfolio

### Product Range

Product Code	Colour	External Catheter Ø		Internal Catheter Ø	Catheter Length	Flow Rate ml/min.
SR+DS2619WD	Violet	26G	0.64 mm	0.45 mm	19 mm	19
SR+DS2419WD	Yellow	24G	0.74 mm	0.55 mm	19 mm	29
SR+DS2225WD	Blue	22G	0.90 mm	0.65 mm	25 mm	42
SR+DS2032WD	Pink	20G	1.00 mm	0.75 mm	32 mm	59
SR+DS1832WD	Green	18G	1.30 mm	0.95 mm	32 mm	103
SR+DS1845WD	Green	18G	1.30 mm	0.95 mm	45 mm	96
SR+DS1745WD	White	17G	1.50 mm	1.15 mm	45 mm	155
SR+DS1645WD	Grey	16G	1.75 mm	1.35 mm	45 mm	225
SR+DS1445WD	Orange	14G	2.00 mm	1.55 mm	45 mm	290



### Packaging

Single Unit Blister pack

Unit Box 50 pieces

Shipping Carton 200 pieces, Dimensions 390 x 384 x 180 mm

### Regulatory Status

CE Certification – Class IIa medical device

#### Notified body

Istituto Superiore di Sanità 0373

#### Manufacturer

Delta Med S.p.A. (Italy)

#### Distributor

TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 Leuven, Belgium

### Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10<sup>-6</sup>

### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper of the unit pack
- Do not store at extreme temperature and humidity

### Precautions

- Do not use if unit package is damaged
- Do not re-use / do not re-sterilize
- Do not attempt to re-insert a partially or completely withdrawn needle
- Dispose of safely after single use to avoid risk of infection
- Application time should be consistent with accepted standards of practice

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.

This range complements the wider Terumo product offering including a comprehensive portfolio of both safety and conventional I.V. products which can be found at our web page on [www.terumo-europe.com](http://www.terumo-europe.com)

Your Terumo partner



#### TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.  
+81 3 6742 8500 +32 16 38 12 11

[www.terumo-europe.com](http://www.terumo-europe.com)

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.