





## mattress replacement system for decubitus prophylaxis and therapy



polyurethane ventilation air cells



Soft-Care-Tex™ cover



also available with bi-stretch cover



- alternating pressure system (19 cm height, about 7.5", incl. 5 cm, 2" foam pocket)
- can be used as a mattress overlay or as a complete system using the insert pocket for a 5 cm (2") foam inlay
- 17 modular air cells made of high-quality polyurethane with 9 laser-perforated cells (ventilation bladders)
- decubitus prophylaxis and therapy up to stage 4 according to EPUAP (European Pressure Ulcer Advisory Panel)
- mattress can be used with an analog or digital pump

### SLK II

max. patient weight: **135 kg (285 lbs)**  
 decubitus stage: **I to IV**  
 low to high risk

### technical specifications

- dimensions mattress: 19 x 90 x 200 cm (7.5" x 35.4" x 78.7")
- dimensions overlay: 14 x 90 x 200 cm (5.5" x 35.4" x 78.7")

### scope of delivery

- SLK II® mattress replacement system with Soft-Care-Tex™ cover (also available with bi-stretch cover) and control unit
- transport bag with carrying straps and zipper, made from durable nylon