

Product Information

VACUETTE® Sodium Fluoride, Potassium Oxalate/ EDTA/ Sodium Heparin; Monoiodoacetate Tubes

Glucose Tubes

Both sodium fluoride and iodoacetate are used in **VACUETTE®** Glucose Tubes as glycolysis inhibitors to preserve glucose when combined with an anticoagulant such as EDTA, potassium oxalate or heparin.

The tubes are used for determinations in the following area: **Glucose testing**

VACUETTE® FE Glucose tubes contain the following per 1 ml blood:

2-4 mg sodium fluoride and 1.2-2 mg K3EDTA

VACUETTE® FH Glucose tubes contain the following per 1 ml blood:

2-4mg sodium fluoride and 12-30 IU sodium heparin

VACUETTE® FX Glucose tubes contain the following per 1 ml blood:

2-4mg sodium fluoride and 1-3mg potassium oxalate

VACUETTE® MJA Glucose tubes contain the following per 1 ml blood:

5.5mg iodoacetate and 12-30 I.U. lithium heparin



Paediatric Tubes

These tubes contain the same additives. The draw volume is either 1 or 2ml of blood and is therefore suitable for children and elderly patients.



VACUETTE® FX Glucose Tubes

Item No.	Nom. Vol.	Cap Colour	Ring Colour	Tube Size	Description
FX Sodium Fluoride / Potassium Oxalate - PREMIUM (Safety Twist Cap)					
454061	2ml	grey	white	13 x 75	FX Sodium Fluoride / Potassium Oxalate, paediatric draw
454094	3ml	grey	black	13 x 75	FX Sodium Fluoride / Potassium Oxalate
454062	4ml	grey	black	13 x 75	FX Sodium Fluoride / Potassium Oxalate
456062	6ml	grey	black	13 x 100	FX Sodium Fluoride / Potassium Oxalate
FX Sodium Fluoride / Potassium Oxalate - Non-ridged (pull cap)					
454238	2ml	grey	white	13 x 75	FX Sodium Fluoride / Potassium Oxalate, paediatric draw
454297	4ml	grey	black	13 x 75	FX Sodium Fluoride / Potassium Oxalate
454438	4ml	grey	black	13 x 75	FX Sodium Fluoride / Potassium Oxalate, transpar. label

Product Information

VACUETTE® Sodium Fluoride, Potassium Oxalate/ EDTA/ Sodium Heparin; Monoiodoacetate Tubes

Glucose Tubes

VACUETTE® FE Glucose Tubes

Item No.	Nom. Vol.	Cap Colour	Ring Colour	Tube Size	Description
FE Sodium Fluoride / K3EDTA - PREMIUM (Safety Twist Cap)					
454085	2ml	grey	white	13 x 75	FE Sodium Fluoride / K3EDTA, paediatric
454091	4ml	grey	black	13 x 75	FE Sodium Fluoride / K3EDTA
FE Sodium Fluoride / K3EDTA - Non-ridged (pull cap)					
454221	2ml	grey	white	13 x 75	FE Sodium Fluoride / K3EDTA, paediatric
454033	4ml	grey	black	13 x 75	FE Sodium Fluoride / K3EDTA

VACUETTE® FH Glucose Tubes

FH Sodium Fluoride / Sodium Heparin - PREMIUM (Safety Twist Cap)					
454218	2ml	grey	white	13 x 75	FH Sodium Fluoride / Sodium Heparin, paediatric
FH Sodium Fluoride / Sodium Heparin - Non-ridged (pull cap)					
454430	2ml	grey	white	13 x 75	FH Sodium Fluoride / Sodium Heparin, paediatric

