



sonicaid®

# Waterproof 2MHz and 3MHz Fetal Dopplers

MkIII

**HUNTLEIGH**  
A MEMBER OF THE ARJO FAMILY

# Sonicaid FD1 & FD3 MkIII

NEW look  
updated range  
of Rate Display  
and Audio  
Dopplers

Fetal Heart Rate Display Dopplers with fixed waterproof 2MHz or 3MHz probe features:

- Compact and light-weight design for ease of portability
- Large high contrast LCD display of FHR in BIG numbers
- Combined Power ON/OFF/Volume control for ease of use
- High sensitivity probes to maximize performance in a wide variety of clinical applications and procedures
- Supports 2MHz and 3MHz waterproof probes suitable for use in water birth
- Easy clean probes for improved infection control
- AA batteries for convenience and long term cost effectiveness



## Technical Specification

Product Features	FD1 MkIII	FD3 MkIII
Product Order Code	FD1 MkIII	FD3 MkIII
Display	FHR Rate Display	FHR Rate Display
Probe Frequency	2MHz	3MHz
Attached Waterproof Probes	•	•
Noise Reduction	Active	Active
Integrated Loudspeaker	•	•
Headphone Output	•	•
Battery Type Supplied as Standard	2 x AA Alkaline batteries	2 x AA Alkaline batteries
Battery Life (no. of 1min exams)	500	500
Warranty*	3 year (1 year on probes, accessories and wires)	
Weight	304gms (10.7oz)	
Dimensions (Main unit)	Height 140mm (5.5") Width 75mm (3.0") Depth 33mm (1.18")	

As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

**Manufactured and distributed by Huntleigh Healthcare Ltd.**  
35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom  
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk  
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:  
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF  
©Huntleigh Healthcare Limited 2019

**A Member of the Arjo Family**

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW: 1001043-3

