


## tergumed 710

### Lateral Flexion



#### TECHNICAL DATA

- » Dimensions (L/W/H): 126/92-115/163 cm
- » Total weight block: 84 kg (4 x 1 kg, 16 x 5 kg)
- » Total weight: up to 370 kg
- » User weight: max. 180 kg
- » Energy consumption standby: 0,0037 kWh/0,0115 kWh\*
- » Energy consumption in operation: 0,018 kWh/0,0266 kWh\*/0,064 kWh\*\*
- » \* With Autolift option
- » \*\* Option Autolift is active

#### SPECIFICATION

- » Strength testing and training for trunk lateral flexion
- » Tested safety in accordance with the European Medical Device Regulation (MDR)
- » smart assist for training control
- » Isometric measurement system
- » Visual feedback strength training with individual curves
- » Automatic seat height adjustment
- » Ultrasonic sensors for weight detection
- » 0 to 84 kg training weight with 1 kg step increments
- » Reproducible scaled positioning
- » Laser pointer for optimal positioning
- » Start position adjustment in 10° steps
- » Stepless trunk fixation adjustment
- » Stepless thigh fixation adjustment
- » Weight block cover for protection

#### FRAME

- » RAL 7021 Black grey frosted

#### CASING

- » RAL 9006 White aluminium

#### UPHOLSTERY

- » F6461197 Anthrazit
- » F6461454 Feuer