

DELIVERY SPECIFICATIONS**Product:** Slabbetje 2L ronde hals, beige 38x70cm**Article:** 2112270

Pcs/crt	600	Pcs/pallet	39600
Crt/row	6	Rows/pallet	11
Crt/pallet	66	Pallet	120x80 cm, One-way
Pcs/bag	100	Bags/crt	6
Dimensions			
	Length [mm]	Width [mm]	Height [mm]
Product	700	380	
Product folded	195	380	
Bag	385	280	
Carton Inner	390	390	200
Carton Outer	398	398	208
Pallet	1200	800	2438
Weights			
	Net	Gross	
Product [g]		9,83	
Bag [g]		9,8	
Bag filled [g]	983,1	992,9	
Carton [g]		450	
Carton filled [kg]	5,90	6,41	
Pallet [kg]	389,3	437,9	
Volume			
Carton filled	32948 cm ³		
Pallet	2,34 m ³		

Note: Product weights may vary up to +/- 15% due to material variation**Material:** Tissue 18 g/m²
Pe 0,014 mm
Reinforcement: Pe 0,014 mm

The tissue is made from virgin pulp


The Pe film is made from polyethylene and meets the conditions to BGA and FDA.

PE-Bag: Neutral white

Carton: Label: Printed as prescript by Customer


Lot number: According to "Lot numbers list"

Product Specification:



Packaging: 600 pcs of folded bibs packed as 100 pcs/bag and 6 bags/cart.

Box:



Pallet



