

Description:	Forceps
Final dimensions:	13 cm
Composition:	Reinforced Polypropylene.
Latex:	Free
PVC:	Free
Re use:	No
Colour:	green
Sterilized by:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.



Nr and Revision Date:	11 -- 2019/06/14
CE Class:	I Steril
CE mark:	CE0120
Directive:	MDD 93/42/EEC amended 2007/47/EC
Applicable norms:	EN ISO 13485:2016 ; EN 1041+A1:2013 ; EN ISO 11607-1:2009 ; EN ISO 11607-2:2006 ; EN ISO 11737-1:2006 ; EN ISO 11737-2:2009 ; EN 556-1:2001 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; EN 1422:2014 ; ISO 1133-1:2011 ; EN 1041:2008 ; ENISO11737-1:2006/AC:2009 ; EN 556-1:2001/AC:2006 ; EN ISO 11135:2014 ; ISO 16142-1:2016 ; EN ISO 15 223-1:2016 ;

Level	Type	Dimension	Stamp	Labels
1	1 Flexible peel-open pack	80x190mm	stamped	
2	35 Folding box	78x193x360mm		labeled
3	560 Carton box	780x380x315mm		labeled

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply
Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
Intended use	Intended to grasp and hold dressing material and other items during medical procedures.

WASTE AND ENVIRONMENT REGULATIONS	Weight (gr) of packing material used in 100 pcs	637
--	---	-----

ISSUED BY: Quality Management Department, 2019/08/08