

# Severo Electric Pill Grinder v3.0

Key Features	Products/Article codes
<ul style="list-style-type: none"> <li>• Fast and effective crushing of medication</li> <li>• Avoiding exposure to medication and medication dust</li> <li>• Small and portable design</li> <li>• Easy to Clean</li> <li>• Near silent operation</li> <li>• Rechargeable battery</li> <li>• Produced in the Netherlands</li> </ul>	Severo Medicine grinder (Art. code: SEVERO3.0), including: <ul style="list-style-type: none"> <li>• Stainless steel grinding head</li> <li>• Recharger</li> </ul> Disposables: <ul style="list-style-type: none"> <li>• Severo cups (Art. code: MED001): 2100 cups per box</li> <li>• Severo rondels (Art. code: MED003): 2100 cups per box</li> <li>• Optional:               <ul style="list-style-type: none"> <li>• Severo cup lids (Art. code: MED002): 2000 lids per box</li> </ul> </li> </ul>

Grinding Specifications	
Grinding type	Autonomous grinding cycles via combined compression and rotation.
Pill grinding	Grinding of multiple standard pills and coated pills per cycle
Cup volume	30 ml
Grinding cycle time	8 seconds for regular pills
Sound level	52 dB (soft whisper level)
Battery last	500 crushing's (= ca three days of continuous use)
Charging time	2 hours
Cost per grinding cycle	0,05 € (when grinding with cups and rondels)
Switches	On/Off
Indicators	Battery load and battery quality indicator
Warranty	2 years on the device, 1 year on the battery (not incl. user error)

Quality Assurance
Before shipment each device is thoroughly tested/inspected on:
1. Scratches, dents, dirt or other exterior deformations
2. Odd sounds during operation
3. Functioning and grinding quality
4. Battery capacity
5. Availability of accessories

Technical Specifications	
Nominal Voltage	12 Volts
Approximate weight	ca 1,5 Kilograms
Dimensions (LxBxH)	80 x 70 x 217 mm
Terminals	5,5 x 2.1 x 12 mm DC Connector for recharging
AC/DC power supply	13 volts at 3 A
Battery	Rechargeable sealed lithium-ion, 12 Volts 3 Ah
Battery self-discharge	ca 1% of capacity per month (25 c°)
Fuse	5 A

