

Product Identification

Product family	Surgical, dental and podiatry instruments
Reference	590062S
Name	Probe 14 cm with graduation 2 mm satin
Sterilization	Ethylene oxide
Brand	Sterisets®



Material

Code	Designation	Qty	Composition
590062	Probe 14 cm with graduation 2 mm satin	1	AISI-410 Stainless steel

Packaging

Type	Sterile single pack
Package weight	4g
Package dimensions	170 x 70 x 10mm
Quantity per inner box	25 pieces
Quantity transport carton	100 pieces

Hazards Identification

Latex	No
Medicinal substances	No
Animal tissue	No
Toxic or dangerous materials	No

Shelf life

5 years from the date of manufacture, unless any component has a shorter shelf-life.

Storage conditions

Keep away from sunlight and keep dry.

Environmental and disposal considerations

Not regarded as dangerous for the environmental. Dispose of material according to local, regional and national regulations.

Compliance Regulations & Standards

Regulation (EU) 2017/745 (as applicable), Directive 93/42/EEC, ISO 13485, ISO 14971, ISO 15223-1, ISO 10993-1, ISO 11737, ISO 11135.

