



Technical data sheet

Product category: Sterile gauze balls

REF: 56203S5

Description: Ball gauze 17 thr. 20 x 17 cm Ø 35 mm with elas 5 pcs/pch *S



Technical drawing:



Material:

Code	Designation	Qty	Specification
56203	Ball gauze 17 thr. 20 x 17 cm Ø 35 mm with elastic ring	5	Absorbent gauze 100% cotton

Package
Weight: 5 grams
Dimensions: 140 x 90 x 10 mm

Sterilization: 121°C Steam

Brand: Sterisets®

Packaging: Sterile single pack

Quantity in pcs	
Inner box	Transport carton
30	120