


Thorax Catheter, bent right-angled



REF: 02.000.10.016, 02.000.10.020, 02.000.10.024,
02.000.10.028, 02.000.10.032, 02.000.10.036

Product group: Sterile Thorax Catheters and Trocar Catheter

Intended use:

The thoracic catheters are inserted into the pleural space via a trocar during surgery to drain fluids and/or air from the intrapleural space.

Product description:

The right-angled thoracic catheters have lateral eyes, a distal opening, and a smooth connector at the proximal end. They also feature an X-ray contrast strip and a 1 cm graduation along the thoracic catheter.

X-ray contrastable: X-ray contrast strip

Scale value: 1 cm

Sterile: Yes

REF	Outer diameter	Length	Number (eyes)	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
02.000.10.016	16 Ch	53 cm	6	10 pcs	04034342074586	100 pcs	04034342262419
02.000.10.020	20 Ch	53 cm	6	10 pcs	04034342263164	100 pcs	04034342263171
02.000.10.024	24 Ch	53 cm	6	10 pcs	04034342074623	100 pcs	04034342263096
02.000.10.028	28 Ch	53 cm	6	10 pcs	04034342074647	100 pcs	04034342263102
02.000.10.032	32 Ch	53 cm	6	10 pcs	04034342074661	100 pcs	04034342262426
02.000.10.036	36 Ch	53 cm	6	10 pcs	04034342074685	100 pcs	04034342263119

Materials

Product	PVC, barium sulphate
Packaging	PE-bag, Tyvek®, PET/PE, cardboard, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

Additional information

Regulatory basis	Directive 93/42/EEC
Notified Body	TÜV SÜD Product Service GmbH (0123)
Classification	Ila
UMDNS	11-308
eClass	34290206
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	EO gas
Shelf life	5 Years
Instructions for Use	Yes