



FM 150 Pro Vein Trainer


 german
 engineering


Air pressure massage of femur and lower legs with inflating and deflating air cushion



Continuously adjustable massage intensity



Improves the venous blood circulation



Revitalising pressure massage

Promotes blood circulation and relieves heavy and tired legs

Looses tension

Air pressure massage with inflating and deflating air cushion

Continuously adjustable massage intensity

With 4 cuffs: 2 femur and 2 lower leg cuffs

Alternate inflation of the leg cuffs

Incl. Safety STOP button

With deflation function - for emptying the air cushions quickly

Incl. timer function 10/20/30 min.

Easy operation with manual control switch

Versatile application - whether you're at home or on the go

Individual adjustment of the leg cuff using hook-and-loop fasteners

Storage pouch included

3 Watt

Medical device

Battery or mains adaptor driven

Mains adaptor voltage: Input: 100-240V~, 50/60Hz, 200 mA
 Output: DC 6V === 500 mA

Incl. 4 x 1.5V battery AA

Product measurements (LxWxH): Femur cuff: 85 x 23,5 x 0,5 cm
 Lower leg cuff: 61 x 67 x 0,5 cm

Product weight: approx. 1,7 kg

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64327 0

Item no.: 643.27