

iHealth Neo

Smart Arm blood pressure monitor BP5S

Compatible with iPhone and Android smartphones



Blood pressure monitoring on iPhone and Android smartphones

Requires iOS v9.0 and + or Android version 4.4 and +
See the list of compatible devices on www.ihealthlabs.eu



iHealth Neo features

- Clinically validated blood pressure monitor
- Systolic, diastolic and heart rate (number of beats per minute)
- 200 measurements memory on the meter
- Irregular Heartbeat Detection (arrhythmia)
- LED display screen
- Tubeless device
- Inflation measurement
- Usable with or without connection
- Save all your measurements in your MyVitals app account



Technical specifications:

- Dimensions: 141,5mm×60,5mm×19mm
- Weight: 240g
- Cuff circumference: 22cm - 42 cm
- Measuring method: Oscillometric method, automatic inflation and measurement
- Connectivity: Bluetooth 4.1 class 2
- Measuring range:
Systolic: 60-260mmHg
Diastolic: 40-199mmHg
Pulse rate: 40-180 beats/minute
- Power supply: 1x 3.7V Li-ion 950mAh
- Storage conditions:
Temp: -20°C~55°C, humidity: ≤90%RH
- Operating Environment:
Temp: 5°C~40°C, humidity: ≤85%RH



In the box:

- 1 iHealth Neo smart blood pressure monitor
- 1 User manual
- 1 Quick Start Guide
- 1 USB cable
- 1 Travel bag

Track your vitals information with the free iHealth MyVitals App



iHealth proposes many other great healthcare products:
activity trackers, scales, blood pressure, SpO2 and blood glucose monitors

CE0197



www.iHealthLabs.eu

iHealth
Connected to my Health