


iHealth ThermoPro

Smart non-contact infrared thermometer NT13B

Compatible with iPhone and Android smartphones



Temperature monitoring on iPhone and Android smartphones

Requires iOS v5.0 and + or Android version 4.3 and +

iHealth ThermoPro features

- Clinically validated
- Easy to use
- Forehead temperature measurement
- Fast reading (1s)
- Automatic switch off
- Bluetooth smart



Technical specifications:

- Connectivity: Bluetooth Smart
Frequency: 2402~2480MHz
Output power range: $\leq 4\text{dBm}$
- Dimensions: 141mm×42mm×55.5mm
- Weight: 120g
- Battery: 2 AAA batteries 1.5V
- Measuring site: Forehead
- Display resolution: 0.1°C
- Range:
Human Body: 34°C~43°C (93.2°F~109.4°F)
Object: 0°C~100°C (32.0°F~212.0°F)
- Accuracy:
Human Body: $\pm 0.2^\circ\text{C} (\pm 0.4^\circ\text{F})$: from 34 to 43 °C (93.2 to 109.4 °F)
Object: $< 40^\circ\text{C} \pm 2^\circ\text{C}; \geq 40^\circ\text{C} \pm 5\%$
- Operating Environment:
Temp: 15°C ~ 35°C, humidity: $\leq 95\%$ RH,
Pressure: 700-1060hPa
- Conditions for transport and storage:
Temp: -25°C ~ 55°C, humidity $\leq 95\%$ RH,
Pressure: 700-1060hPa

iHealth

Compatible with iHealth SDK

iHealth proposes many other great healthcare products: blood pressure monitors, blood glucose meters, pulse oximeter, scales and activity trackers.



In the box:

- 1 iHealth ThermoPro thermometer
- 1 User guide
- 2 AAA batteries

CE0197

www.iHealthLabs.eu

iHealth[™]
Connected to my health