

HOSPITAL CARE

AGANI™

Hypodermic needle
Product Information



Sterile hypodermic needle for single use.

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect¹.

Description	Raw Materials
1 Hub	Polypropylene + master batch
2 Cannula	Stainless steel
Adhesive	Epoxy glue
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.

Sterilization
With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the packaging
- Do not store at extreme temperature and humidity

Precautions

- Check integrity of the package before use
- Do not reuse
- Do not use if unit package is damaged
- Dispose of safely after single use to avoid risk of infection

Your needles of choice for a wide range of applications

Product Range

Product Code	Colour	Cannula ø	Cannula Length		Bevel	
AN*1838R1	Pink	18G	1.20 mm	1 ½ inch	38 mm	Regular 11°
AN*1838S1	Pink	18G	1.20 mm	1 ½ inch	38 mm	Short 17°
AN*1925R1	Cream	19G	1.10 mm	1 inch	25 mm	Regular 11°
AN*1938R1	Cream	19G	1.10 mm	1 ½ inch	38 mm	Regular 11°
AN*1938S1	Cream	19G	1.10 mm	1 ½ inch	38 mm	Short 17°
AN*1950R1	Cream	19G	1.10 mm	2 inch	50 mm	Regular 11°
AN*2025R1	Yellow	20G	0.90 mm	1 inch	25 mm	Regular 11°
AN*2038R1	Yellow	20G	0.90 mm	1 ½ inch	38 mm	Regular 11°
AN*2116R1	Green	21G	0.80 mm	5/8 inch	16 mm	Regular 11°
AN*2125R1	Green	21G	0.80 mm	1 inch	25 mm	Regular 11°
AN*2138R1	Green	21G	0.80 mm	1 ½ inch	38 mm	Regular 11°
AN*2150R1	Green	21G	0.80 mm	2 inch	50 mm	Regular 11°
AN*2232R1	Black	22G	0.70 mm	1 ¼ inch	32 mm	Regular 11°
AN*2238R1	Black	22G	0.70 mm	1 ½ inch	38 mm	Regular 11°
AN*2316R1	Blue	23G	0.60 mm	5/8 inch	16 mm	Regular 11°
AN*2325R1	Blue	23G	0.60 mm	1 inch	25 mm	Regular 11°
AN*2332R1	Blue	23G	0.60 mm	1 ¼ inch	32 mm	Regular 11°
AN*2425R1	Purple	24G	0.55 mm	1 inch	25 mm	Regular 11°
AN*2516R1	Orange	25G	0.50 mm	5/8 inch	16 mm	Regular 11°
AN*2525R1	Orange	25G	0.50 mm	1 inch	25 mm	Regular 11°
AN*2613R1	Brown	26G	0.45 mm	½ inch	13 mm	Regular 11°
AN*2623R1	Brown	26G	0.45 mm	7/8 inch	23 mm	Regular 11°
AN*2716R1	Grey	27G	0.40 mm	5/8 inch	16 mm	Regular 11°
AN*2719R1	Grey	27G	0.40 mm	¾ inch	19 mm	Regular 11°
AN*3013R1	Yellow	30G	0.30 mm	½ inch	13 mm	Regular 11°



Packaging

Single unit
Blister pack (Film/Paper)

Unit box 100 pieces
Dimensions 105 x 84 x 86 mm
For item numbers AN*1950R1 and AN*2150R1
Dimensions 105 x 84 x 110 mm

Shipping carton 10.000 pieces
Dimensions 450 x 440 x 450 mm
For item numbers AN*1950R1 and AN*2150R1
5.000 pieces,
Dimensions 545 x 240 x 450 mm

Regulatory Status

CE Certification –
Class IIa medical device

Notified body TÜV Product Service GmbH, TÜV SUD Group 0123

Manufacturer Zhejiang Kindly Medical Devices Co. Ltd. (China)

EC representative Shanghai International Corp. GmbH (Europe) - Hamburg, Germany

Distributor TERUMO EUROPE N.V. Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a wide range from 18G up to 30G

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

Your Terumo partner



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.