


iHealth Gluco+

Smart glucose monitor BG5s

Compatible with iPhone and Android smartphones



Glucose monitoring on iPhone and Android smartphones

Requires iOS v9.0 and + or Android version 4.4 and +
See the list of compatible devices on www.ihealthlabs.eu

iHealth Gluco+ features

- Measurement and monitoring of glucose levels for diabetic patients
- Instantly displays glucose levels on the meter or smartphone
- 500 measurements memory on the meter
- Bright LED screen for easy reading
- Tracks data, graphs and history on the iHealth Gluco-Smart App
- Bluetooth 4.1 BLE sync, easy installation
- Compatible only with iHealth test strips
- Evaluation of measurements as recommended by international guidelines
- Share your medical data with you healthcare professional or family
- Charges by mini USB
- Secure and free data storage



Technical specifications:

- Dimensions: 98mmx35mmx27.8mm
- Weight: 32.5g
- Measuring method: Amperometric technology using glucose dehydrogenase (FAD-GDH)
- Measuring range: 20~600 mg/dL, 1.1 ~ 33.3mmol/L
- Power supply: DC 3.7V Li-Ion 250mAh
- Storage conditions:
 - Test Strips: Temp: 4°C~30°C, humidity <85%RH
 - Blood glucose monitor: Temp: -20°C~55°C, humidity <85%RH
- Operating Environment: Temp: 10°C~40°C
- Blood Source: fresh whole capillary blood
- Volume of blood: min 0.5µl
- For self-measurement only



In the box:

- 1 iHealth Gluco+ glucose monitor
- 1 quick start guide
- 1 user manual
- 1 USB charging cable
- 1 travel case



Track your vitals information with the free iHealth Gluco-Smart App



iHealth proposes many other great healthcare products:
Activity trackers, scales, blood pressure monitors and pulse oximeter

CE0197

www.iHealthLabs.eu

iHealth
Connected to my Health