



CETRO® Hook

9601




CETRO MEDICAL


CETRO® Hook är en hinnsprängare avsedd för amniotomi i samband med förlossning. Den är utformad för att vara enkel att använda och ge ett stabilt grepp vid handhavande. CETRO® Hook har en rundad och skyddad spets i kombination med en vass hake för att möjliggöra en kontrollerad öppning av fosterhinnorna. Längden är anpassad för god åtkomst vid användning.

• Greppvänlig design • Rundad och skyddad spets • Vass hake • Anpassad längd	• Steril • Latexfri
---	--

MD CE MDM Plastic Lebanon, Bldg 1682, Al Salam Village Exit, Al Shamaa Blvd., Haret Saidia., Saidia - Lebanon.



ARTIKEL NR.	PRODUKT	AVD. FÖRPACKNING	TRANSP.FÖRPACKNING	TRANSP.PACKNING
9601	7,6 g	100 st	1600 st	13 kg



CETRO® Hook

CETRO® Hook is an amniotic membrane perforator intended for use during amniotomy in labor. It is designed to be easy to use and provide a stable grip during handling. CETRO® Hook features a rounded and protected tip combined with a sharp hook to enable controlled perforation of the amniotic membrane. The length is adapted to allow good access during use.

The product is supplied sterile and is intended for single use.

- Easy-grip design
 - Rounded and protected tip
 - Sharp hook
 - Adapted length
-
- Sterile
 - Latex Free

Material: Polystyrene impact resistant (opaque)

MD

CE












STERILE R



MDM Plastic Lebanon, Bldg 1682, Al Salam Village Exit, Al Shamaa Blvd., Haret Saidaa,, Saida - Lebanon.



ARTICLE NO.	PRODUCT	DEPT.PACKAGING	TRANSP.PACKAGING	TRANSP.PACKAGING
9601	7,6 g	100 pcs	1600 pcs	13 kg

