




Garment Measuring Guide

The following guide is intended as an aid to help the user select the correct size garment.

The treatment type will determine the garment type e.g. hydroven 1 for uniform compression or hydroven 3 for graduated pressure.

When determining the correct size garment:



^{*}Measurements should be made at the widest point on the calf and thigh.

Hydroven 3 Arm

Three Chamber Arm Garments



Order Code	Туре	Length (L)	Shoulder Circumference (S)	Hand Circumference (H)
1533360	Full Arm	68cm	62cm	44cm
1533361	Full Arm	78cm	62cm	44cm

Hydroven 3 Leg

Three Chamber Leg Garments



Order Code	Туре	Length (L)	Thigh Circumference (T)	Calf Circumference (C)
1533333	Half Leg	50cm	61cm	56cm
1533334	Full Leg	66cm	64cm	56cm
1533335	Full Leg	71cm	66cm	56cm
1533336	Full Leg	76cm	72cm	56cm
1533337	Full Leg	84cm	72cm	56cm
1533338	Full Leg	92cm	72cm	56cm

Hydroven Inserts

Compatible with Hydroven 1 and Hydroven 3 garments



Order Code	Туре	Length (L)	Added Circumference Wide End	Added Circumference Narrow End
		Leg Inserts		
1533362	Half Leg Insert	50cm	19cm	14cm
1533363	Full Leg Insert	66cm	19cm	14cm
1533364	Full Leg Insert	71cm	19cm	14cm
1533365	Full Leg Insert	76cm	19cm	14cm
1533366	Full Leg Insert	84cm	19cm	14cm
1533367	Full Leg Insert	92cm	19cm	14cm
Arm Inserts				
1533368	Full Arm Insert	68cm	17cm	12cm
1533369	Full Arm Insert	78cm	17cm	12cm

Hydroven 1 Arm

Single Chamber Arm Garments



Order Code	Туре	Length (L)	Shoulder Circumference (S)	Hand Circumference (H)
1533317	Half Arm	51cm	56cm	44cm
1533331	Full Arm	68cm	62cm	44cm
1533332	Full Arm	78cm	62cm	44cm

Hydroven 1 Leg

Single Chamber Leg Garments



Order Code	Туре	Length (L)	Thigh Circumference (T)	Calf Circumference (C)
1533311	Half Leg	50cm	61cm	50cm
1533312	Full Leg	66cm	64cm	54cm
1533313	Full Leg	71cm	66cm	55cm
1533314	Full Leg	76cm	72cm	56cm
1533315	Full Leg	84cm	72cm	57cm
1533316	Full Leg	92cm	72cm	57cm

As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excelle nee and improved performance, for life.

Manufactured and distributed by H untleigh Healthcare Ltd. 35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office: ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfords hire, LU5 SXF @Huntleigh Healthcare Limited 2019

A Member of the Arjo Family @ and " are trademarks of Huntleigh Technology Limited As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW:: 1001043-2



ARSEUS MEDICAL

Official distributor Huntleigh

Arseus Medical

Rijksweg 10 - 2880 Bornem

0800 76 769 - info.hospital@arseus-medical.be

www.arseus-medical.be

