



HK 25

beurer
wellbeing

HK 25 Heat pad

 german
engineering



Textile cover machine washable
up to 40°C



Electronic temperature control

3 illuminated temperature settings

Automatic switch off after approx. 90 minutes

BSS (Beurer Safety System)

Detachable control

Oeko-Tex

The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.

Mains operation:

Input: 220-240V~; 50-60 Hz; 100 W

Product size (LxB): 40 x 30 cm

Product weight: approx. 460 g

3 year guarantee

Sales unit: 10 / Shipping carton: -

EAN no.: 4211125 25502 2

Item no.: 255.02