

- Device classification: Class 1
- Complies with the 93/42/CEE directive and its amendments
- GMDN code: 40811



FIAB SpA - Via P. Costoli, 4 - 50039 VICCHIO - FIRENZE - ITALY

<i>General description</i>	Adapter for snap electrodes	
<i>Attachments</i>	Yes	Photo
<i>Batch reference</i>	No	
<i>Connected device</i>	Yes	Electrocardiograph
<i>Single use device</i>	No	
<i>Sterile device</i>	No	
<i>Sterilisation</i>	No	
<i>Autoclavable</i>	No	
<i>Cleaning</i>	Yes	Can be washed with a dry cloth dipped in alcohol.
<i>Unsuitable substance</i>		
<i>Component parts</i>	No	
<i>Material</i>		
<i>Labelling</i>	Yes	COMED, reference, name, address, EC, website
<i>Packaging</i>	Yes	Individually in a plastic bag.
<i>Special packaging</i>	Yes	The temperature should be between 0°C and + 50°C and relative humidity between 20% and 80%.
<i>Duration of use</i>	Average	
<i>Marking on medical device</i>	No	
<i>References</i>	26 700 24	
<i>Designation</i>	Adapter for snap electrodes	
<i>Dimensions L x W x H (mm)+- 3mm</i>	35 x 18 x 10	
<i>Net weight (g)+1%</i>	5	
<i>Colours</i>	White	
<i>Connection</i>	For snap fastener	
<i>Temperature resistance (°C)</i>	0°C to +50°C	
<i>Usage</i>	Device for connecting standard cables different types to different electrodes for medical use (ECG, electrotherapy...)	
<i>Instruction for use</i>	Connect the adapter to the electrode by touch.	
<i>Recycling</i>	No	
<i>Recovery</i>	In medical waste containers	
<i>WEEE/REACH</i>	No/Yes	
<i>Phthalates/Latex</i>	No/No	
<i>RoHS</i>	Yes	
<i>Safety data</i>	R-Phrases = None Combination phrases = None	

COMED is not responsible marketing and does not have manufacturer status according to Directive 93/42/EEC, acting solely as a distributor.