



## Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of  
ENRAF-NONIUS B.V., Vareseweg 127, 3047 AT Rotterdam, The Netherlands

Quality Management System according to ISO 13485:2016, EN ISO 13485:2016  
and Annex VII of European Medical Devices Directive 93/42/EEC.

Device name:	Manumed Series (NPU)
Description:	Treatment couches for physical and manual therapy
Article number(s):	See appendix CFG-5500.000-48
Classification:	I (according to rule 1, Annex IX of MDD 93/42/EEC)
Record of conformity:	026-400.269-4x
Device marking:	CE

We hereby declare that the above-mentioned device  
complies with the European Medical Devices Directive 93/42/EEC and it has been notified to the Dutch  
Ministry of Health.

Any modifications to the product not authorized by Enraf-Nonius will invalidate this declaration.

Rotterdam, date: 28 February 2019

Signature:

A handwritten signature in blue ink, appearing to be 'Th.J. Doodkorte', written over a horizontal line.

Th.J. Doodkorte  
Commercial Director



Configuration code of Manumed  
Appendix CFG-5500000-48  
Page 1 of 2

**N212N.N004.L5R.00.00**

type                      frame colour   padding                      upholstery colour                      width                      inlorface   ew                      voltage                      packing

**Type**

L001 Basic-2 section. head  
L002 Basic. 2-section. back  
L004 Traction. Fixed height  
L005 Bobath. 1-section. Fixed Height  
N2xx Optimal 3-section  
N3xx Optimal 2-section  
N4xx Osteo  
N5xx Bobath 1- & 2- section  
N6xx Comfort 3-section  
N7xx Special tilt  
N9xx Traction

**Frame colour**

G Chrome  
N Grey  
P White  
U Black

**Padding**

N normal  
H extra hard  
S extra thick and soft  
L Luxury

**Upholstery colours**

smooth skai  
012 Black  
019 Bordeaux  
127 Yellow  
135 Beige  
138 Mint  
139 Signal Red  
141 Turquoise  
155 Jade  
169 Cornflower Blue  
199 Lavender  
204 Sirius Blue  
213 Grey

**extra rough skai**

811 Black  
812 Grey  
851 Beige  
852 Red  
888 White  
933 Blue

**Section Width**

L5R 57 cm rectangular nose hole. straight leg  
L6R 67 cm rectangular nose hole. straight leg  
L6T 67 cm rectangular nose hole. trapezium leg  
L8T 80 cm rectangular nose hole. trapezium leg  
D5R 57 cm with triangle nose hole. straight leg  
D6R 67 cm with triangle nose hole. straight leg  
D6T 67 cm with triangle nose hole. trapezium leg  
D8T 80 cm with triangle nose hole. trapezium leg  
N6R 67 cm with small nose hole. straight leg  
Z0R 100 cm without nose hole. straight leg  
Z2R 120 cm without nose hole. straight leg  
Z5R 57 cm without nose hole. straight leg  
Z6R 67 cm without nose hole. straight leg  
Z6T 67 cm without nose hole. trapezium leg  
Z8T 80 cm without nose hole. trapezium leg

Rotterdam, date: 28 February 2019

Signature:

Th.J. Doodkorte  
Commercial Director



Configuration code of Manumed  
Appendix CFG-5500000-48  
Page 2 of 2

**Interface**

- 0 hydraulic
- A foot switch
- B 2 pc foot switch
- C hand switch
- K table switch
- L 4 button hand switch
- M 4 button table switch
- N foot switch + table switch
- R all-around (bar ↓ - couch ↓)
- S all around + table switch (bar ↓ - couch ↓)
- T all-around (bar ↑ - couch ↓)
- U all-around + table switch (bar ↑ - couch ↓)

**ex works options**

- 0 no options
- 1 fixation rails
- 2 prepared hole in seat
- 3 ergo wings
- 4 foot / wheel extension 30 mm
- 5 options 4+1
- 6 options 4+2
- 7 options 4+3
- 8 options 4+1+2
- 9 options 1+2

**Mains voltage**

- C 230V CH 8.000N motor
- D 230V CH 10.000N motor
- E 230V EUR 8.000N motor
- F 230V EUR 10.000N motor
- G 230V UK 8.000N motor
- H 230V UK 10.000N motor
- J 115V 8.000N motor
- K 115V 10.000N motor
- 0 hydraulic pump or fixed height

**Packaging**

- 0 cardboard top cover
- 1 box pallet
- 2 cardboard top cover + access
- 3 carton box + carton blocks
- 4 fumigated pallet + box
- D cardboard box

Rotterdam, date: 28 February 2019

Signature:



Th.J. Doodkorte  
Commercial Director