

EU Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of **Enraf-Nonius B.V.**
Quality Management System according to EN ISO 13485:2016

Product name:	Manumed
Catalogue number(s):	5500000
	Delivery according to Manumed configurator 103-5500000-4A on page 2.
Basic UDI-DI:	871999250855000003A
Intended purpose:	The Manumed is a treatment couch, which supports a patient during treatment by a physiotherapist
Risk class:	I, according to rule 13 in Annex VIII of (EU) 2017/745
Device marking:	CE

The device covered by this declaration, in combination with its test certificate is in
conformity with:

(EU) 2017/745 (Medical Device Regulation)

2006/42/EC (Machinery Directive)

2011/65/EU as amended by (EU) 2015/863 (RoHS Directive)

Any modifications to the product not authorized by Enraf-Nonius will invalidate this
declaration.

Signed for and on behalf of Enraf-Nonius B.V.:

T. Doodkorte
Commercial Director

Place: Rotterdam
Date of issue: 30 April 2021

Enraf-Nonius B.V.
Vareseweg 127
3047 AT Rotterdam
The Netherlands
SRN (NL-MF-000000429)

EU Declaration of Conformity

Configuration code of Manumed

Appendix 103-5500000-4A

N212N.N004.L5R.00.00

type frame colour padding upholstery colour width interface ew voltage packing

Type

L001	Basic-2 section, head
L002	Basic, 2-section, back
L004	Traction, Fixed height
L005	Bobath, 1-section, Fixed Height
N2xx	Optimal 3-section
N3xx	Optimal 2-section
N4xx	Osteo
N5xx	Bobath 1- & 2- section
N6xx	Comfort 3-section
N7xx	Special tilt
N9xx	Traction

Frame colour

G	Chrome
N	Grey
P	White
U	Black

Padding

N	normal
H	extra hard
S	extra thick and soft
L	Luxury

Upholstery colours

012	Black
019	Bordeaux
127	Yellow
135	Beige
138	Mint
139	Signal Red
141	Turquoise
155	Jade
169	Cornflower Blue
170	Baby Blue
199	Lavender
204	Sirius Blue
213	Grey
406	Petrol
407	Cherry
408	Terracotta
409	Skin
410	Anthracite
811	Volcano Black
812	Grey
850	Walnut
851	Sand
852	Red
886	Leaf
888	Snow
933	Blue

Section Width

L5R	57 cm rectangular nose hole, straight leg
L6R	67 cm rectangular nose hole, straight leg
L6T	67 cm rectangular nose hole, trapezium leg
L8T	80 cm rectangular nose hole, trapezium leg
D5R	57 cm with triangle nose hole, straight leg
D6R	67 cm with triangle nose hole, straight leg
D6T	67 cm with triangle nose hole, trapezium leg
D8T	80 cm with triangle nose hole, trapezium leg
N6R	67 cm with small nose hole, straight leg
Z0R	100 cm without nose hole, straight leg
Z2R	120 cm without nose hole, straight leg
Z5R	57 cm without nose hole, straight leg
Z6R	67 cm without nose hole, straight leg
Z6T	67 cm without nose hole, trapezium leg
Z8T	80 cm without nose hole, trapezium leg

Interface

0	hydraulic
A	foot switch
B	2 pc foot switch
C	hand switch
K	table switch
L	4 button hand switch
M	4 button table switch
N	foot switch + table switch
R	all-around (bar ↓ - couch ↓)
S	all around + table switch (bar ↓ - couch ↓)
T	all-around (bar ↑ - couch ↓)
U	all-around + table switch (bar ↑ - couch ↓)

Mains voltage

C	230V CH 8.000N motor
D	230V CH 10.000N motor
E	230V EUR 8.000N motor
F	230V EUR 10.000N motor
G	230V UK 8.000N motor
H	230V UK 10.000N motor
J	115V 8.000N motor
K	115V 10.000N motor
O	hydraulic pump or fixed height
V	230V EUR + safety switch
W	230V EUR 10.000N + safety switch

ex works options

0	no options
1	fixation rails
2	prepared hole in seat
3	ergo wings
4	foot / wheel extension 30 mm
5	options 4+1
6	options 4+2
7	options 4+3
8	options 4+1+2
9	options 1+2

Packaging

0	cardboard top cover
1	box pallet
2	cardboard top cover + access
3	carton box + carton blocks
4	fumigated pallet + box
D	cardboard box

Enraf-Nonius B.V.

Vareseweg 127

3047 AT Rotterdam

The Netherlands

SRN (NL-MF-000000429)