

# TECHNICAL DATA SHEET



## Cotton balls



MATERIAL	<ul style="list-style-type: none"> <li>• 100 % bleached absorbent cotton or cotton/viscose mix</li> <li>• H2O2 bleached (non chlorine bleached, free of any optical brightener)</li> <li>• Meeting the European Pharmacopoeia Standards</li> <li>• Oeke-Tex 100 certified</li> <li>• In accordance with REACH European Standards</li> <li>• Bioburden analysis on request</li> </ul>
PRODUCT INFORMATION	<ul style="list-style-type: none"> <li>• Used for absorption of exudates, peripheral cleaning wounds and supporting the application of products on the skin + laboratory plugs EO (Ethylene Oxide) sterilizable</li> <li>• Different sizes available: from 0,3 to 5 gr</li> <li>• Shelf life: 5 years from the date of manufacture</li> <li>• Class I medical device, single use product</li> <li>• General tolerance: +/- 10 %</li> </ul>
PACKAGING	<ul style="list-style-type: none"> <li>• Packed in polyethylene bags to protect from dust and humidity</li> <li>• The goods can also be delivered in bulk</li> <li>• Other packaging and packaging available on request</li> <li>• * Available in Cotoneige brand</li> <li>• Export box sizes: 60x40x60 cm (16) and 59x39x24 cm (36)</li> <li>• Should be stored in a clean and dry place, away from the sunlight</li> </ul>

### CONDITIONNEMENTS STANDARDS

Reference	Description	Balls/bag	Weight/Bag	Bags/Box	Box/Pallet
1500B03	1500 balls 0.3 gr	1500	(450 gr)	12	16
100B06COT*	100 balls 0.6 gr	100	(60 gr)	25	36
100BO06COTC*	100 coloured balls 0.6 gr	100	(60 gr)	25	36
500B06	500 balls 0,6-0,7gr	500	(325 gr)	16	16
700B06	700 balls 0,6 gr – 0,7 gr	700	(455 gr)	16	16
500G03	500g of balls 0,3 gr	(1666)	500 gr	10	16
500G06	500g of balls 0,6 gr	(833)	500 gr	10	16
500G09	500g of balls 0,9 gr	(555)	500 gr	10	16
5KG06	5 kg of balls 0,6-0,7 gr	(7692)	5000 gr	1	16
5KG2	5 kg of balls 2 gr	(2500)	5000 gr	1	16
5KG3	5 kg of balls 3 gr	(1666)	5000 gr	1	16
5KG5	5 kg of balls 5 gr	(1000)	5000 gr	1	16

#### ALAN & CO

Avenue André Ernst 17 – 4800 Verviers – Belgium  
Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69

We are ISO 13485 certified and our products are CE marked

Version 2021/02