

PRODUCT DATA SHEET

ProSet Perican®(4457927)

Set for Epidural Anesthesia

Attributen Naam	Attribuut Waarde
Product line	RA
Type of set	Complete set
Intended use	Single shot epidural anaesth.
Bestseller	No
Epidural Needle Bevel	Tuohy
Epidural Needle Gauge	18
Epidural length	80.00 mm
Catheter-Type	-
Catheter-Gauge	-
Catheter-Eyes	-
Perifix Filter	-
Catheter connector	-
PinPad	-
Catheter fixation	No

Attributen Naam	Attribuut Waarde
Safety Components	No
Syringes: LS or LL	Luer-Slip
incl. LOR Syringe	Perifix LOR 8 ml LS
incl. Syringes	1 x Omnifix 3 ml - LS, 1 x Omnifix 5 ml - LS, 1 x Omnifix 10 ml - LS
incl. Needles	1 x Sterican 21G 0.8x40 mm, 1 x Sterican 18G 1.2x40 mm, 1 x Sterican 26G 0.45x25 mm
Plasters	No
Type of Plaster	-
includes drugs	-
Gauze compresses	Yes
Gauze compresses: Size	5x Gauze co. 7.5x7.5 cm white
Swabs	No
Swabs: Size	-
Swab forceps	No
Type of Swab forceps	-
Scissors	-
Swabsticks	Yes
Swabsticks: Quantity	2
Patient drape	No
Patient drape: Size	-
Patient drape: Type	-
Patient drape: Color	-
Patient drape: Window	No
Drape: Window size	-
Patient drape: Opening size	-
Patient drape: Adhesive tape	No

Attributen Naam	Attribuut Waarde
Patient drape: Adhesive pad	No
Extra surgical drape	1
Surgical drape: Size	75x90 cm
Surgical drape: Color	blue
Surgical drape: Adhesive tape	Yes
Field sheet: Quantity	1
Field sheet: Size	75x90 cm
Field sheet	Yes
Ultrasound cover	No
Ultrasound cover: Size	-
Ultrasound gel	No
Extra Disinfection tray	1x 265 ml white
Protective clothing	No
Protective clothing: Type	-
Suture material	No
Suture material: Type	-
Thread Material	-
Size of Thread (EP)	-
Angle of Needle 1	-
Drug preparation	No
Drug preparation: Type	-
incl. scalpel	No
Sterility	sterile
Product contains Latex	No
Disposable / Reusable	disposable

Attributen Naam	Attribuut Waarde
NRFit	No