

WIEGAND[®] PillCrusher Professional



Art. Nr. 202.95.01
PillCrusher Professional



Art. Nr. 202.95.11
Disposable Pouches for
PillCrusher Professional,
50 pcs.

With the WIEGAND[®] PillCrusher Professional, tablets can be processed quickly, hygienically and without applying force into a ready-to-use powder.



WIEGAND[®] PillCrusher Professional

- No contamination of the environment, or personnel by drug residues or dust.
- Depending on the type of pills, simultaneously several pills can be crushed together.
- Protects joints due to minimal effort.
- Ideal for tube feeding, because tablets are crushed into a fine powder which doesn't clog the tubes.
- The WIEGAND[®] PillCrusher Professional is silent in use and doesn't disturb the patients.
- The waveform in the pressing area facilitates the crushing of large or hard tablets.
- Note: Use only with WIEGAND[®] PillCrusher Pouches

Disposable PillCrusher Pouches

- The strong ribbed reinforced PillCrusher Pouches prevent tearing, despite heavy pressures while crushing.
- The pouches seams are bevelled to ensure complete emptying (reduction of the dosage by loss of substance).

5 years Warranty

The WIEGAND[®] PillCrusher Professional is made of steel and impact-resistant plastic. He is extremely durable, easy to clean and maintenance free. Wiegand offers a 5 year warranty on the PillCrusher Professional.

Please consider the instructions before use.