



# Stand scale

7830

- ▶ Comfortable display height
- ▶ Solid platform surface
- ▶ Instant-start-function by stepping on the scale
- ▶ Easy to clean surface

## Features

---

The personal scale with stand is the ideal scale for fitness studios, pharmacies, care or other social institutions. When stepping on the scale, the display activates automatically and the weight is shown with red backlit digits. The tread surface is large and stable, but can still be easily transported with the optional set of castors. The scale can be operated via mains power or optionally with a rechargeable battery and can thus also be used flexibly.

- ▶ Bright luminous LED display, 15 mm digit height
- ▶ Extra large tread: 320 x 400 mm
- ▶ Instant-start-function by stepping on the scale
- ▶ Tare function
- ▶ Mother-child-weighing
- ▶ Easy to clean surface
- ▶ No assembly required - the scale is immediately ready for use
- ▶ Service interface (resp. data cable available)
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Mains operation
- ▶ Optional rechargeable battery operation
- ▶ Optional wheels for the easy transport of the scale to the patient

## Products

---

Order number Variation	Approvable	Division g	Maximum load kg	Platform size mm
7830.30.001 Stand scale	No	100	250	365 x 490

---

## Accessories

---

Order number Variation	Description
2550.03.019	Data cable for service interface
2557.12.002	Rechargeable battery, installed (must be included in the purchase order, no retrofit possible)
2513.12.002	Wheels (2 pieces), installed (must be included in the purchase order, no retrofit possible)