

HOSPITAL CARE

AGANI™

Safety hypodermic needle

Product Information



The Agani™ Safety hypodermic needles, when fitted with a standard luer slip or luer lock syringe are intended for the aspiration and injection of fluids for medical purposes. After withdrawing the needle from the body, the attached needle safety shield can be activated to cover the needle immediately after use to minimize the risk of an accidental needle stick.

	Description	Raw materials
1	Sheath	Polypropylene
2	Collar	Polypropylene
3	Hub	Polypropylene
4	Cannula	Stainless steel
	Protector	Polypropylene (PP)
	Adhesive	Epoxy glue
	Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.



Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper side of the blister package
- Do not store at extreme temperature and humidity

Secure safety activation

Reliable sharps protection with one-hand activation

Product range

Product Code	Colour	Cannula ø	Cannula Length	Bevel
SAN1825R1	Pink	18G 1.2 mm	1 inch 25 mm	Regular 11°
SAN1838R1	Pink	18G 1.2 mm	1 ½ inch 38 mm	Regular 11°
SAN1925R1	Cream	19G 1.1 mm	1 inch 25 mm	Regular 11°
SAN1938R1	Cream	19G 1.1 mm	1 ½ inch 38 mm	Regular 11°
SAN2025R1	Yellow	20G 0.9 mm	1 inch 25 mm	Regular 11°
SAN2038R1	Yellow	20G 0.9 mm	1 ½ inch 38 mm	Regular 11°
SAN2125R1	Green	21G 0.8 mm	1 inch 25 mm	Regular 11°
SAN2138R1	Green	21G 0.8 mm	1 ½ inch 38 mm	Regular 11°
SAN2225R1	Black	22G 0.7 mm	1 inch 25 mm	Regular 11°
SAN2238R1	Black	22G 0.7 mm	1 ½ inch 38 mm	Regular 11°
SAN2325R1	Blue	23G 0.6 mm	1 inch 25 mm	Regular 11°
SAN2338R1	Blue	23G 0.6 mm	1 ½ inch 38 mm	Regular 11°
SAN2516R1	Orange	25G 0.5 mm	5/8 inch 16 mm	Regular 11°
SAN2525R1	Orange	25G 0.5 mm	1 inch 25 mm	Regular 11°
SAN2538R1	Orange	25G 0.5 mm	1 ½ inch 38 mm	Regular 11°
SAN2613R1	Brown	26G 0.45 mm	½ inch 13 mm	Regular 11°
SAN2713R1	Grey	27G 0.4 mm	½ inch 13 mm	Regular 11°
SAN3013R1	Yellow	30G 0.3 mm	½ inch 13 mm	Regular 11°



Not all item codes will be available in phase 1 of launch.
Please contact your local Terumo contact for more information.

Choose from a wide range from 18G up to 30G

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

Packaging

Single unit Blister pack (film/paper)

Unit box 100 pieces
Dimensions 215 x 195 x 85 mm

Shipping carton 800 pieces
Dimensions 455 x 415 x 205 mm

Precautions

- Do not use if package is damaged
- Do not reuse
- Use once and discard immediately in accordance with local safety standards

Regulatory Status

CE Certification –
Class IIa medical device

Notified body
TÜV Product Service GmbH,
TÜV SÜD Group 0123

Manufacturer
Zhejiang Kindly Medical Devices
Co. Ltd. (China)

EC representative
Shanghai International Corp. GmbH
(Europe) - Hamburg, Germany

Distributor
TERUMO EUROPE N.V.,
Interleuvenlaan 40, 3001 Leuven,
Belgium

Your Terumo partner



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.