

tergumed 710

Rotation

TECHNICAL DATA

- » Dimensions (L/W/H): 128-130/148/163-171 cm
- » Total weight block: 84 kg (4 x 1 kg, 16 x 5 kg)
- » Total weight: up to 400 kg
- » User weight: max. 180 kg
- » Energy consumption standby: 0,0037 kWh/0,0115 kWh*
- » Energy consumption in operation: 0,018 kWh/0,0266 kWh*/0,064 kWh**
- » * With Autolift option
- » ** Option Autolift is active



SPECIFICATION

- » Strength testing and training for trunk rotation
- » Tested safety in accordance with the European Medical Device Regulation (MDR)
- » smart assist for training control
- » Isometric measurement system
- » Visual feedback strength training with individual curves
- » Automatic seat height adjustment
- » Ultrasonic sensors for weight detection
- » 0 to 84 kg training weight with 1 kg step increments
- » Reproducible scaled positioning
- » Start position adjustment in 10° steps
- » Stepless shoulder fixation adjustment
- » Stepless knee fixation adjustment
- » Stepless footrest adjustment
- » Stepless hip fixation adjustment
- » Weight block cover for protection

FRAME

- » RAL 7021 Black grey frosted

CASING

- » RAL 9006 White aluminium

UPHOLSTERY

- » F6461197 Anthrazit
- » F6461454 Feuer