

HOSPITAL CARE

Obturator

for I.V. Catheter

Product Information



To be used in
combination with I.V. catheters

Serves to close the catheter and to prevent clot formation during intermittent infusion.

An obturator with the same colour code as the catheter is available for each size.

Description	Raw Materials
Obturator	Polypropylene
Protector	Polyethylene

Does not contain components made of PVC or natural rubber latex.



Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10^{-6}



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper of the blister pack
- Do not store at extreme temperature and humidity



Precautions

- Do not use if unit pack is damaged
- Do not re-use / do not resterilize
- Dispose of safely after single use to avoid risk of infection

Easy closure for all types of I.V. catheters

Product Range

Product Code	Colour	Obturator Size		Length
SO+DM2419WD	Yellow 	24G	0.74 mm	19 mm
SO+DM2419PD	Yellow 	24G	0.74 mm	19 mm
SO+DM2225WD	Blue 	22G	0.90 mm	25 mm
SO+DM2225PD	Blue 	22G	0.90 mm	25 mm
SO+DM2032WD	Pink 	20G	1.00 mm	32 mm
SO+DM2032PD	Pink 	20G	1.00 mm	32 mm
SO+DM1832WD	Green 	18G	1.30 mm	32 mm
SO+DM1832PD	Green 	18G	1.30 mm	32 mm
SO+DM1845WD	Green 	18G	1.30 mm	45 mm
SO+DM1845PD	Green 	18G	1.30 mm	45 mm
SO+DM1745WD	White 	17G	1.50 mm	45 mm
SO+DM1745PD	White 	17G	1.50 mm	45 mm
SO+DM1645WD	Grey 	16G	1.75 mm	45 mm
SO+DM1645PD	Grey 	16G	1.75 mm	45 mm
SO+DM1445WD	Orange 	14G	2.00 mm	45 mm
SO+DM1445PD	Orange 	14G	2.00 mm	45 mm



WD = for I.V. Catheter without injection port:
SURFLO™ S-Plus, VERSATUS S-Plus, SURFLO™-W,
VERSATUS-W, Surshield™ VERSATUS-S/-W
PD = for I.V. Catheter with injection port:
SURFLO™-WP, VERSATUS-WP, Surshield™ VERSATUS-WP

Packaging

Single Unit Blister pack

Unit Box 50 pieces

Shipping Carton

500 pieces

Gross Weight ± 2.2 kg

Dimensions

± 720 x 200 x 190 mm

Regulatory Status

CE Certification –
Class IIa medical device

Notified body

TÜV SÜD Product Service GmbH (0123)

Manufacturer

Delta Med S.p.A. (Italy)

Distributor

TERUMO EUROPE N.V., Interleuvenlaan
40, 3001 Leuven, Belgium

This range complements the wider Terumo product offering including a comprehensive portfolio of both safety and conventional I.V. products which can be found at our web page on www.terumo-europe.com

Your Terumo partner



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

Image for visualisation purpose only