

This site is intended for Healthcare Professionals. Not all products are available in all countries. Please contact your local office for more information. If you are a patient, it is important to discuss health information with your doctor. To continue, please acknowledge and dismiss.

Accept & Dismiss Notification Disclaimer

[Back to Ethicon Site](#) →

[Home](#) / [Product Code Search](#) | Ethicon Product Center

**ETHICON**


*Johnson & Johnson* SURGICAL TECHNOLOGIES

# TVC55

## Product Name

PROXIMATE® Linear Cutters

Quick View	
Product Category	Surgical Stapling
Product Group	Linear Cutter



Additional Specifications	
Description	55mm Vascular Linear Cutter with White Reload
Brand	PROXIMATE®
Staple Technology	B-Form
Device Technology	Cam Closure Mechanism
Size	55mm
Max Firings	8
Color	WHITE
Staple Rows	4
Tissue Thickness	Vascular/Thin

# of Staples	56
Cut Line Length	53mm
Staple Line Length	57mm
Category	Surgical Stapling
QTY/BX	3

## DISCLAIMER

Before using any medical device, review all relevant package inserts with particular attention to the indications, contraindications, warnings and precautions, and steps for use of the device(s). We are committed to providing you with the most recent and accurate product information and will update this site on an as needed basis. The images used on the website are for reference only and does not reflects the actual configuration of the product. We do not provide any guarantee as to the accuracy, completeness or suitability of the images shown of the products for any particular purpose. You acknowledge that such information and materials may contain inaccuracies or errors and we expressly exclude liability for any such inaccuracies or errors to the fullest extent permitted by law.

© Ethicon US, LLC. 2020.  
085133-171129

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

This site is published by **Ethicon US, LLC**, which is solely responsible for its contents, and is intended for U.S. audiences only.

This site is intended for Healthcare Professionals. If you are a patient, it is important that you discuss

[News](#)

[Careers](#)

[CA state requirements](#)

[Contact us](#)

[Johnson & Johnson Institute](#)

[Legal notice/Terms and conditions](#)

[Privacy policy](#)

[Cookie Policy](#)

[Supplemental contract terms](#)

[Choice Matters](#)

[The Microwave Ablation](#)



**Product Catalog**

**Order Product**

information about the benefits and risks of products with your doctor.

Last updated March 26, 2020. The third-party trademarks used herein are trademarks of their respective owners. Data and usage fees may apply. See your carrier for more information.

**Movement**

**Do not sell or share my personal information**

**Customize Cookie Settings**