

HYDROVEN[®] 3 SERIES

3 Chamber Intermittent Pneumatic Compression





An effective treatment for a variety of circulatory disorders

Venous leg ulcers (VLU) are a significant health problem that afflicts a large percentage of the population.

Intermittent Pneumatic Compression (IPC) is an effective treatment for a variety of circulatory disorders¹. The hydroven 3 Intermittent Compression System (IPC) has been developed to address the following vascular conditions:

- Oedema
- Lymphoedema
- Chronic Venous Insufficiency
- Acute and chronic wound management including venous leg ulcers and post-surgical wounds.

hydroven 3 can be used as an adjunct to traditional compression therapies in patients who may be immobile or are unable to tolerate pain from conventional therapy or have previously failed treatment with other modalities. The simple to use system is designed to provide each of its patients with maximum comfort and a treatment which delivers an effective outcome.



"Hydroven 3 is a simple to use, lightweight pump with preset cycle times and adjustable pressure settings for individual comfort"



hydroven 3 can be used at home, clinic or hospital setting, its robust design ensures that the system can be handled, transported and stored safely and easily.

hydroven 3 offers clinicians a choice of single or three chamber garments which can be applied to the patients' upper or lower limbs. Garment inflation and deflation cycles which are specific to the type of therapy being delivered are controlled from a high reliability silent pneumatic pump.

A legacy you can trust

Huntleigh has been innovating and developing products that assist in the management and treatment of Vascular and Lymphatic conditions for over 30 years. Linking this experience to our world class hand held Dopplers and assessment systems, Huntleigh are uniquely positioned to offer clinicians comprehensive solutions for the holistic assessment and treatment of vascular conditions.

Effective treatment with si

Snap lock connection system

A secure single hose, snap-lock connection system ensures a positive connection between both garment and pump. It eliminates any potential for accidental disconnection or air leak.

Silent operation

hydroven 3 uses a vibration free pump which minimises disturbances around the patient environment, ensuring that the system fits comfortably into the patients' lifestyle and surrounding area.

Single limb or bi-lateral assessment

Single or bi-lateral assessment is possible by simply connecting two garments to the system.

mple, easy to use operation

Integrated handle

The integrated handle promotes safe and effective transportation.



Simple controls

The simple minimal controls can be accessed with the minimum of effort. Main on / off control and garment pressure adjustment promotes patient compliance which enhances the patient experience.

Auto detection

Automatic detection of single or dual garments supports the different treatment regimes.

Designed for patient comfort, durability and performance

Durable garments

All hydroven garments are manufactured with a soft pliable inner lining for maximum patient comfort; the construction of each garment ensures inward expansion to fit the contours of each limb. The inflatable foot section ensures pressure is distributed around the foot to aid blood flow and treatment. If patient movement is restricted or impaired, the large ring pull zip ensures easy garment application and removal.



Clinical evidence

Application of Intermittent Pneumatic Compression (IPC) in conjunction with sustained graduated compression improves outcomes in patients with the most advanced venous diseases¹.

Using IPC in combination with bandaging may be more effective than bandaging alone. The sequential inflation and deflation of the chambers creates intermittent

pressure peaks, mimicking the effect of the calf muscle pump thus offering a number of benefits such as improved tissue oxygenation, accelerated venous return and healing rates².

Management of chronic venous leg ulcers is complex and demanding. Application of IPC has shown to be beneficial and could be considered to optimise outcome³.

Specification

hydroven 3 Intermittent Pneumatic Compression System			
Order code	5100XX (XX defined country variant)		
Supply voltage	110VAC to 240VAC		
Supply frequency	50Hz		
Pressure range	Pressure Range - 30 - 100mmHg +/- 5%		
Power input	14VA		
Size	270 x 140 x 140mm (10.6 x 5.5 x 5.5")		
Weight	2.4Kg (5.3lbs)		
Degree of protection against electric shock	Class II, double insulated Type BF		
Degree of protection against liquid ingress	IPXO – no protection		
Mode of operation	Continuous		
Order Code	Description	Size	

	Single Chamber Leg Garments	
5101L50	Half leg garment	500 mm (19.5")
5101L66	Full leg garment	660mm (26")
5101L71	Full leg garment	710mm (28")
5101L76	Full leg garment	760mm (30")
5101L84	Full leg garment	840mm (33")
5101L92	Full leg garment	920mm (36.25")
	Single Chamber Arm Garments	
5101A51	Half arm garment	510 mm (20")
5101A68	Full arm garment	680 mm (26.75")
5101A78	Full arm garment	780 mm (30.6")
	Three Chamber Leg Garments	
5103L50	Half leg garment	500 mm (19.6")
5103L66	Full leg garment	660 mm (26")
5103L71	Full leg garment	710 mm (28")
5103L76	Full leg garment	760 mm (30")
5103L84	Full leg garment	840 mm (33")
5103L92	Full leg garment	920 mm (36.25")
	Three Chamber Arm Garments	
5103A68	Full arm garment	680 mm (26.75")
5103A78	Full arm garment	780 mm (30.6")
	Leg Inserts	
510LI50	Half leg insert	500 mm (19.6")
510LI66	Full leg insert	660 mm (26")
510LI71	Full leg insert	710 mm (28")
510LI76	Full leg insert	760 mm (30")
510LI84	Full leg insert	840 mm (33")
510LI92	Full leg insert	920 mm (36.25")
	Arm inserts	
510AI68	Full arm insert	680 mm (26.75")
510AI78	Full arm insert	780 mm (30.6")

References

1 Comerota, A J. Intermittent Pneumatic Compression: Physiologic and clinical basis to improve management of venous leg ulcers. Journal of Vascular Surgery 53, 4 2011 2 World Union of Wound Healing Societies (WUWHS). Principles of best practice. Compression in venous leg ulcers. A consensus document. London: MEP Ltd, 2008. 3 Case Reports in Clinical Medicine, 2014, 3, 513-517. Published on line September 2014. http://www.scrip.org/journal/crcm

As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

Manufactured and distributed by Huntleigh Healthcare Ltd. 35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office: ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF ©Huntleigh Healthcare Limited 2019

A Member of the Arjo Family

® and ^m are trademarks of Huntleigh Technology Limited As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW:: 1001043-2







We would like to express our sincere appreciation to the Welsh Wound Innovation Centre for their assistance and for allowing access to their facility and resources.

HUNTLEIGH