

EU Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of **Enraf-Nonius B.V.** Quality Management System according to EN ISO 13485:2016

Product name:	ManuMax
Catalogue number(s):	1570000
	Delivery according to ManuMax configurator 103-1570000-44 on page 2.
Basic UDI-DI:	871999250815700003K
Intended purpose:	The ManuMax is a treatment couch, which supports a patient during treatment by a physiotherapist
Risk class:	I, according to rule 13 in Annex VIII of (EU) 2017/745
Device marking:	CE

The device covered by this declaration, in combination with its test certificate is in conformity with:

(EU) 2017/745 (Medical Device Regulation)
2006/42/EC (Machinery Directive)
2011/65/EU as amended by (EU) 2015/863 (RoHS Directive)

Any modifications to the product not authorized by Enraf-Nonius will invalidate this declaration.

Signed for and on behalf of Enraf-Nonius B.V.:

T. Doodkorte
Commercial Director

Place: Rotterdam
Date of issue: 30 April 2021

Enraf-Nonius B.V.
Vareseweg 127
3047 AT Rotterdam
The Netherlands
SRN (NL-MF-000000429)

EU Declaration of Conformity

Configuration code of ManuMax Appendix 103-157000-44

TTT	F	S	C	Z	A	HL	F	P	CCC	WW	N	V	M	P	OOOO
Type	Folds	Sections	Frame colour	Side cover colour	Actuator	High Low	Flexion	Uphol. padding	Uphol. colour	Section width	Breathing hole	Mains voltage	Mains Plug	Packaging	Options
TTT	Type			F	Flexion				139	Signal red				OOOO	Options
EXA	Examination			N	No				019	Wine red				0	No
MUL	Multi			H	Gas spring				852	Red				1	Fixation rail for armpit / shoulder supports
EXE	Exercise			T	Table switches				933	Blue				2	Fixation rail for gyno supports
TIL	Tilt								513	Grey Crib 5				3	Holder for Pelvic support
TRA	Traction			P	Upholstery padding				WW	Section width				4	Paperroll holder foot section
MAN	Manipulation			N	Normal				57	57 cm					
				L	Luxe				67	67 cm					
F	Folds			CCC	Upholstery colour				80	80 cm					
1-4	Number of sections in longitudinal direction			886	Leaf				10	100 cm					
				887	Fuchsia				12	120 cm					
S	Sections			888	Snow				N	Breathing hole					
1-9	Total number of sections			851	Sand				T	Triangular					
				779	Walnut				L	Long					
C	Frame colour			812	Antracite				N	No					
G	Grey			811	Volcano				V	Mains Voltage					
				204	Sirius blue				S	100-240 V					
Z	Side cover colour			169	Cornflower blue				M	Mains Plug					
D	Dark			170	Baby blue				E	Europe					
L	Light			199	Lavender				U	United Kingdom					
				213	Grey				J	Japan					
A	Actuator			012	Black				A	North America					
H	Hydraulic			406	Cherry				C	Switzerland					
E	Electric			407	Ocre				P	Packaging					
				408	Petrol				S	Standard					
HL	High/Low operation			409	Skin										
AA	All Around			410	Terracotta										
FF	2x Foot Switch			127	Yellow										
TS	Table Switches			135	Beige										
HS	Hand Switch			138	Mint										
				155	Jade										
				141	Turquoise										

TTT	F	S	C	Z	A	HL	F	P	CCC	WW	N	V	M	P	OOOO																							
Type	Folds	Sections	Frame colour	Side cover colour	Actuator	HL operation	Flexion	Upholstery padding	Upholstery colour	Section width	Breathing hole	Mains	Plug	Packaging	Options																							
			G Grey	D Dark	L Light	H Hydraulic	E Electric	AA All Around	FF 2x Foot Switch	TS Table Switches	HS Hand Switch	N No	G Gas spring	T Table Switches	N Normal	L Luxe	CCC Upholstery colour	55 55 cm	67 67 cm	80 80 cm	10 100 cm	12 120 cm	T Triangular	L Long	N No	S 100-240 V	E Eur	U UK	J Japan	A North America	C CH / Switzerland	S Standard	0 no options	1 Fixation rail for armpit /	2 Fixation rail for gyno set	3 holder for pelvic support	4 Paperroll holder foot section	
Examination		ManuMax, Examination, two fold, 2 section, back section	S	S	X		S	S	X	X	X	S			S			S	S	X	X	X	X	S	S												X	
Multi Three		ManuMax, Multi, three fold, 3 section, head section	S	S	X		S	S	X	X	X	S			X	S		S	S	X	X	X	X	S	S													X
Multi	3	ManuMax, Multi, three fold, 5 section, head section with adjustable armrests	S	S	X		S	S	X	X	X	S			X	S		S	S	X	X	X	X	S	S													X

S=Standard
X=Optional

Enraf-Nonius B.V.
Vareseweg 127
3047 AT Rotterdam
The Netherlands
SRN (NL-MF-000000429)