



PO 60

**beurer**  
medical

## PO 60 *Bluetooth*<sup>®</sup> Pulse Oximeter

 **german  
engineering**

 **Bluetooth**<sup>®</sup>  
**beurer**  
CONNECT 



Optimum value control (SpO<sub>2</sub> and pulse) via app



Determining arterial oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse frequency)

Innovative networking between smartphone and pulse oximeter. Measurements can be monitored using Beurer HealthManager app/Internet/PC

Transfer via *Bluetooth*<sup>®</sup> technology

Incl. 100 memory spaces

Simple and completely pain-free measurement

Small and light – for use at home and on the move

Particularly suitable for persons with:

- Heart failure
- Chronic obstructive pulmonary diseases
- Bronchial asthma

Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)

Easy-to-read colour display

Colour display with 4 display formats

Adjustable display brightness

The display orients automatically (vertical format, horizontal format)

Graphical display of pulse

Low battery indicator

Automatic switch-off

Incl. retaining strap and belt pouch

Medical product

Incl. 2 x AAA batteries

Product measurements (LxWxH): 58.5 x 30 x 32 mm

Product weight: approx. 49 g (including batteries)

5 year guarantee


Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4211125 45420 3


Item no.: 454.20



optional



## CONNECTION



Innovative networking between smartphone and pulse oximeter - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via *Bluetooth®*.



Syncs with  
 • iOS ≥ 10.0  
 • Android™ ≥ 5.0  
 • Bluetooth® ≥ 4.0



List of supported devices

\* tested online application „beurer HealthManager“ Cloud  
 Can only be used with its full functionality with the „beurer HealthManager“ app.

454.20\_PO60\_2019-01-25\_01\_DS\_EN\_BEU Subject to errors and changes

The *Bluetooth®* word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.  
 Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.  
 Google Play and the Google Play logo are trademarks of Google LLC.  
 Android is a trademark of Google LLC.