

compass 540

Leg Extension/Curl



TECHNICAL DATA

- » Dimensions (L/W/H): 117/116-160/167 cm
- » Total weight block: 80 kg (16 x 5 kg)
- » Total weight: up to 270 kg
- » User weight: max. 250 kg
- » Energy consumption standby: 0,0037 kWh
- » Energy consumption in operation: 0,018 kWh

SPECIFICATION

- » Strength training for legs
- » Tested safety in accordance with the European Medical Device Regulation (MDR)
- » 0 to 80 kg training weight with 5 kg step increments
- » Reproducible scaled positioning
- » Start position adjustment from the training position
- » Stepless lever arm length adjustment
- » Stepless seat depth adjustment
- » Knee fixation adjustment
- » Stepless thigh support adjustment
- » Darts frost cover with surface for storage
- » Training panel with smart operating information on the device

OPTIONS

- » smart assist
- » Isometric measurement system
- » Additional weight increments 4 x 1 kg
- » Range of motion limiter
- » 100 kg weight block (20 x 5 kg)
- » Lumbar cushion

FRAME

- » Engine silver
- » RAL 9003 Signal-white (optional)
- » RAL 7021 Black grey frosted (optional)

UPHOLSTERY

- » F6461197 Anthrazit
- » F6461454 Feuer

COVERING

- » darts frost