

HOSPITAL CARE

Obturator

for I.V. Catheter

Product Information



To be used in combination with I.V. catheters

Serves to close the catheter and to prevent clot formation during intermittent infusion.

An obturator with the same colour code as the catheter is available for each size.

| Description | Raw Materials |
|-------------|---------------|
| Obturator | Polypropylene |
| Protector | Polyethylene |

Does not contain components made of PVC or natural rubber latex.



Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10^{-6}



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper of the blister pack
- Do not store at extreme temperature and humidity











Precautions

- Do not use if unit pack is damaged
- Do not re-use / do not resterilize
- Dispose of safely after single use to avoid risk of infection

Easy closure for all types of I.V. catheters

Product Range

| Product Code | Colour | Obturator Size | Length |
|--------------|--|----------------|--------|
| SO+DM2419W | Yellow  | 24G 0.74 mm | 19 mm |
| SO+DM2419P | Yellow  | 24G 0.74 mm | 19 mm |
| SO+DM2225W | Blue  | 22G 0.90 mm | 25 mm |
| SO+DM2225P | Blue  | 22G 0.90 mm | 25 mm |
| SO+DM2032W | Pink  | 20G 1.00 mm | 32 mm |
| SO+DM2032P | Pink  | 20G 1.00 mm | 32 mm |
| SO+DM1832W | Green  | 18G 1.30 mm | 32 mm |
| SO+DM1832P | Green  | 18G 1.30 mm | 32 mm |
| SO+DM1845W | Green  | 18G 1.30 mm | 45 mm |
| SO+DM1845P | Green  | 18G 1.30 mm | 45 mm |
| SO+DM1745W | White  | 17G 1.50 mm | 45 mm |
| SO+DM1745P | White  | 17G 1.50 mm | 45 mm |
| SO+DM1645W | Grey  | 16G 1.75 mm | 45 mm |
| SO+DM1645P | Grey  | 16G 1.75 mm | 45 mm |
| SO+DM1445W | Orange  | 14G 2.00 mm | 45 mm |
| SO+DM1445P | Orange  | 14G 2.00 mm | 45 mm |



W = for I.V. Catheter without injection port:
 SURFLO™ S-Plus, VERSATUS™ S-Plus, SURFLO™-W,
 VERSATUS™-W, Surshield™ VERSATUS™-S/-W
 P = for I.V. Catheter with injection port:
 SURFLO™-WP, VERSATUS™-WP, Surshield™ VERSATUS™-WP

Packaging

Single Unit Blister pack

Unit Box 50 pieces

Shipping Carton

500 pieces

Gross Weight 4.6 kg

Dimensions

720 x 200 x 190 mm

Regulatory Status

CE Certification –
 Class IIa medical device

Notified Body

TÜV Rheinland 0197

Manufacturer

TERUMO EUROPE N.V.

Interleuvenlaan 40, 3001 Leuven, Belgium

This range complements the wider Terumo product offering including a comprehensive portfolio of both safety and conventional I.V. products which can be found at our web page on www.terumo-europe.com

Your Terumo partner



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
 +81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.