

compass 600

Leg Extension/Curl



TECHNICAL DATA

- » Dimensions (L/W/H): 113/132-148/165 cm
- » Total weight block: 80 kg (16 x 5 kg)
- » Total weight: up to 412 kg
- » User weight: max. 200 kg
- » Energy consumption standby: 0,0034 kWh
- » Energy consumption in operation: 0,018 kWh/0,022 kWh*
- » * Fully equipped (position and pulse sensors)

SPECIFICATION

- » Strength training for legs
- » Tested safety in accordance with the European Medical Device Regulation (MDR)
- » smart assist for training control
- » Integrated RFID or chip card system
- » Visual feedback strength training with individual curves
- » Ultrasonic sensors for weight detection
- » 0 to 84 kg training weight with 1 kg step increments
- » Reproducible scaled positioning
- » Start position adjustment from the training position
- » Stepless lever arm length adjustment
- » Stepless seat depth adjustment
- » Knee fixation adjustment
- » Laser pointer for optimal positioning
- » Weight block cover for protection

OPTIONS

- » Isometric measurement system
- » Position and pulse sensor
- » Range of motion limiter
- » Laser pointer for optimal positioning
- » 100 kg weight block (20 x 5 kg)
- » Lordosis pad

FRAME

- » RAL 7021 Black grey frosted

CASING

- » RAL 9006 White aluminium

UPHOLSTERY

- » F6461197 Anthrazit
- » F6461454 Feuer