



PO 80

beurer
medical

PO 80 Pulse Oximeter


 **german**
engineering



Continuous recording of values (SpO₂ and pulse) for up to 24 hours



PC software for detailed display and evaluation



SOFTWARE



Determining arterial oxygen saturation (SpO₂) and heart rate (pulse frequency)

Simple and completely pain-free measurement

Small and light – for use at home and on the move

Particularly suitable for persons with:

- Heart failure
- Chronic obstructive pulmonary diseases
- Bronchial asthma
- Sleep apnoea

Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)

Continuous recording of up to 24h

PC software "beurer SpO₂ Assistant" for detailed display and evaluation

Operating system: Windows 8.1 or higher

USB cable and USB charger

Can be charged by USB (integrated lithium-ion battery)

Adjustable display brightness

Graphical display of pulse and optional pulse tone

Alarm function – limit values can be individually set

Battery status indicator

Automatic switch-off

Incl. retaining strap, belt pouch and metal storage box

Medical device

Product measurements: 57 x 32 x 30 mm

Product weight: approx. 42 g

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 421112545440 1

Item no.: 454.40