



## tergumed 710 Cervical Trainer

### TECHNICAL DATA

- » Dimensions (L/W/H): 142/83-126/196 cm
- » Total training weight: 1-18 kg (0,1 kg steps)
- » Total weight: up to 320 kg
- » User weight: max. 180 kg
- » Energy consumption standby: 0,0037 kWh/0,0115 kWh\*
- » Energy consumption in operation: 0,018 kWh/0,0266 kWh\*/0,064 kWh\*\*
- » \* With Autolift option
- » \*\* Option Autolift is active

### SPECIFICATION

- » Strength testing and training for cervical extension/flexion and lateral flexion
- » Tested safety in accordance with the European Medical Device Regulation (MDR)
- » smart assist for training control
- » Isometric measurement system
- » Visual feedback strength training with individual curves
- » Automatic seat height adjustment
- » Ultrasonic sensors for weight detection
- » 1 to 18 kg training weight with 0,1 kg step increments
- » Reproducible scaled positioning
- » Laser pointer for optimal positioning
- » Start position adjustment in 7,5° steps
- » Stepless lever arm length adjustment
- » Stepless trunk fixation adjustment
- » Stepless thigh support adjustment
- » Weight block cover for protection

### FRAME

- » RAL 7021 Black grey frosted

### CASING

- » RAL 9006 White aluminium

### UPHOLSTERY

- » F6461197 Anthrazit
- » F6461454 Feuer