

## Trocar Catheter



**REF:** 02.000.30.008, 02.000.30.010, 02.000.30.012, 02.000.30.016, 02.000.30.020, 02.000.30.024, 02.000.30.028, 02.000.30.032

**Product group:** Sterile Thorax Catheters and Trocar Catheter

**Intended use:**

The thoracic catheters are inserted into the pleural space via a trocar to drain fluids and/or air from the intrapleural space.

**Product description:**

The thoracic catheters have an internal trocar as an insertion aid, lateral eyes, and a distal opening. They also feature an X-ray contrast strip and a 1 cm graduation along the thoracic catheter.

Sizes Ch 8, 10, 12, 16 have a funnel attachment at the proximal end.

Sizes Ch 20, 24, 28, 32 have a smooth connector.

**X-ray contrastable:** X-ray contrast strip

**Sterile:** Yes

REF	Outer diameter	Length	Number (eyes)	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
02.000.30.008	8 Ch	24 cm	2	10 pcs	04034342365080	100 pcs	04034342365097
02.000.30.010	10 Ch	24 cm	2	10 pcs	04034342074869	100 pcs	04034342263133
02.000.30.012	12 Ch	24 cm	2	10 pcs	04034342074883	100 pcs	04034342263140
02.000.30.016	16 Ch	26 cm	3	10 pcs	04034342074906	100 pcs	04034342348120
02.000.30.020	20 Ch	40 cm	3	10 pcs	04034342074920	100 pcs	04034342262495
02.000.30.024	24 Ch	40 cm	3	10 pcs	04034342074944	100 pcs	04034342262501
02.000.30.028	28 Ch	40 cm	3	10 pcs	04034342074968	100 pcs	04034342262518
02.000.30.032	32 Ch	40 cm	3	10 pcs	04034342074982	100 pcs	04034342273903

## Materials

Product	PVC, PP, barium sulphate, trocar: aluminium (from size Ch 20), stainless steel (Ch 8-16)
Packaging	PE-bag, Tyvek®, PET/PE, cardboard, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

## Additional information

Regulatory basis	Directive 93/42/EEC
Notified Body	TÜV SÜD Product Service GmbH (0123)
Classification	Ila
UMDNS	11-308
eClass	34290206
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	EO gas
Shelf life	5 Years
Instructions for Use	Yes