

Silicone balloon catheter type Foley

REF:

23.000.14.208 - Ch 8
 23.000.14.210 - Ch 10
 23.000.14.212 - Ch 12
 23.000.14.214 - Ch 14
 23.000.14.216 - Ch 16
 23.000.14.218 - Ch 18
 23.000.14.220 - Ch 20

23.000.14.222 - Ch 22
 23.000.14.224 - Ch 24
 23.000.24.216 - Ch 16
 23.000.24.218 - Ch 18
 23.000.24.220 - Ch 20
 23.000.24.222 - Ch 22
 23.000.24.224 - Ch 24

Unit: 10 pcs.

Packing Unit: 500 pcs.

Description

The catheters are used for controlled bladder release and post-operative drainage.

REF	CH	length	volume of balloon	colour code	guide needle
23.000.14.208	8	30 cm	3 ml	black	yes
23.000.14.210	10	30 cm	3 ml	grey	yes
23.000.14.212	12	40 cm	5-10 ml	white	no
23.000.14.214	14	40 cm	5-10 ml	green	no
23.000.14.216	16	40 cm	5-10 ml	orange	no
23.000.14.218	18	40 cm	5-10 ml	red	no
23.000.14.220	20	40 cm	5-10 ml	yellow	no
23.000.14.222	22	40 cm	5-10 ml	purpule	no
23.000.14.224	24	40 cm	5-10 ml	blue	no
23.000.24.216	16	40 cm	5-10 ml	orange	no
23.000.24.218	18	40 cm	20 ml	red	no
23.000.24.220	20	40 cm	30 ml	yellow	no
23.000.24.222	22	40 cm	30 ml	purple	no
23.000.24.224	24	40 cm	30 ml	blue	no

23.000.14.208-224; 23.000.24.212-224

QM/Rev.-Nr.00/26.06.2003

QM/Rev.-Nr.01/25.02.2005

QM/Rev.-Nr.02/09.07.2008

QM/Rev.-Nr.03/16.06.2009

QM/Rev.-Nr.04/14.03.2013

QM/Rev.-Nr.05/27.01.2014

QM/Rev.-Nr.06/21.10.2015

QM/Rev.-Nr.07/31.05.2016

QM/Rev.-Nr.08/02.06.2016

QM/Rev.-Nr.09/07.06.2016

QM/Rev.-Nr.10/24.03.2017

QM/Rev.-Nr.11/13.04.2017

QM/Rev.-Nr.12/05.08.2019

QM/Rev.-Nr.13/27.09.2019

Material Composition

<u>Product:</u>	silicone, ABS, barium sulfate
<u>guide needle:</u>	PBT
<u>Packaging:</u>	individual packing:
	box:
	outer case:
	PE bag, paper, PET/PE
	cardboard
	corrugated cardboard

Sterilisation

EO-Gas

23.000.14.208-224; 23.000.24.212-224

QM/Rev.-Nr.00/26.06.2003

QM/Rev.-Nr.01/25.02.2005

QM/Rev.-Nr.02/09.07.2008

QM/Rev.-Nr.03/16.06.2009

QM/Rev.-Nr.04/14.03.2013

QM/Rev.-Nr.05/27.01.2014

QM/Rev.-Nr.06/21.10.2015

QM/Rev.-Nr.07/31.05.2016

QM/Rev.-Nr.08/02.06.2016

QM/Rev.-Nr.09/07.06.2016

QM/Rev.-Nr.10/24.03.2017

QM/Rev.-Nr.11/13.04.2017

QM/Rev.-Nr.12/05.08.2019

QM/Rev.-Nr.13/27.09.2019