

# Technical drawings

## Commode/shower chair soft seat, self-propelled


Article no:

LI2138.011 - Commode/shower chair soft seat, self-propelled with closed front, 610 mm

LI2138.111 - Commode/shower chair soft seat, self-propelled with closed front, 530 mm

LI2138.021 - Commode/shower chair soft seat, self-propelled with open front, 610 mm

LI2138.121 - Commode/shower chair soft seat, self-propelled with open front, 530 mm



	A	B	C	D
LI2138.011/LI2138.021	570 mm/22.4"	495 mm/19.5"	550 mm/21.7"	320-480 mm/12.6"-18.9"
LI2138.111/LI2138.121	490 mm/19.3"	415 mm/16.3"	470 mm/18.5"	320-400 mm/12.6"-15.7"

	mm	inch
E	1080	42.5"
F	450	17.7"
G	415	16.3"
H	490	19.3"
I	250	9.8"
J	650	25.6"
K	450	17.7"